

PDHA2 siRNA (h): sc-89026

BACKGROUND

The pyruvate dehydrogenase (PDH) complex is a nuclear-encoded mitochondrial matrix enzyme complex that functions as the primary link between glycolysis and the tricarboxylic acid (TCA) cycle by catalyzing the irreversible conversion of pyruvate into acetyl-CoA, an essential step in aerobic glucose metabolism. PDHA2 (pyruvate dehydrogenase α 2), also known as PDHAL, is a 388 amino acid mitochondrial matrix protein expressed in postmeiotic spermatogenic cells. Composed of a tetramer containing two α and two β subunits, PDHA2 consists multiple copies of three enzymatic components: pyruvate dehydrogenase (E1), dihydrolipoamide acetyltransferase (E2) and lipoamide dehydrogenase (E3). PDHA2 is suggested to participate in cell proliferation and may be involved in prostate cancer. PDHA2 is encoded by a gene located on human chromosome 4q22.3, which encodes nearly 6% of the human genome and has the largest gene deserts (regions of the genome with no protein encoding genes) of all of the human chromosomes.

REFERENCES

- Dahl, H.H., Brown, R.M., Hutchison, W.M., Maragos, C. and Brown, G.K. 1990. A testis-specific form of the human pyruvate dehydrogenase E1 α subunit is coded for by an intronless gene on chromosome 4. *Genomics* 8: 225-232.
- Brown, R.M., Dahl, H.H. and Brown, G.K. 1990. Pyruvate dehydrogenase E1 α subunit genes in the mouse: mapping and comparison with human homologs. *Somat. Cell Mol. Genet.* 16: 487-492.
- Online Mendelian Inheritance in Man, OMIM™. 1990. Johns Hopkins University, Baltimore, MD. MIM Number:179061. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Fitzgerald, J., Hutchison, W.M. and Dahl, H.H. 1992. Isolation and characterization of the mouse pyruvate dehydrogenase E1 α genes. *Biochim. Biophys. Acta* 1131: 83-90.
- Lanneluc, I., Mulsant, P., Saidi-Mehtar, N. and Elsen, J.M. 1996. Synteny conservation between parts of human chromosome 4q and bovine and ovine chromosomes 6. *Cytogenet. Cell Genet.* 72: 212-214.
- Klar, J., Asling, B., Carlsson, B., Ulvsbäck, M., Dellsen, A., Ström, C., Rhedin, M., Forslund, A., Anneren, G., Ludvigsson, J.F. and Dahl, N. 2005. RAR-related orphan receptor A isoform 1 (ROR α 1) is disrupted by a balanced translocation t(4;15)(q22.3;q21.3) associated with severe obesity. *Eur. J. Hum. Genet.* 13: 928-934.
- Kumar, V., Rangaraj, N. and Shivaji, S. 2006. Activity of pyruvate dehydrogenase A (PDHA) in hamster spermatozoa correlates positively with hyperactivation and is associated with sperm capacitation. *Biol. Reprod.* 75: 767-777.
- Kumar, V., Kota, V. and Shivaji, S. 2008. Hamster sperm capacitation: role of pyruvate dehydrogenase A and dihydrolipoamide dehydrogenase. *Biol. Reprod.* 79: 190-199.
- Pinheiro, A., Faustino, I., Silva, M.J., Silva, J., Sá, R., Sousa, M., Barros, A., de Almeida, I.T. and Rivera, I. 2010. Human testis-specific PDHA2 gene: methylation status of a CpG island in the open reading frame correlates with transcriptional activity. *Mol. Genet. Metab.* 99: 425-430.

CHROMOSOMAL LOCATION

Genetic locus: PDHA2 (human) mapping to 4q22.3.

PRODUCT

PDHA2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PDHA2 shRNA Plasmid (h): sc-89026-SH and PDHA2 shRNA (h) Lentiviral Particles: sc-89026-V as alternate gene silencing products.

For independent verification of PDHA2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-89026A, sc-89026B and sc-89026C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

PDHA2 siRNA (h) is recommended for the inhibition of PDHA2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PDHA2 gene expression knockdown using RT-PCR Primer: PDHA2 (h)-PR: sc-89026-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.