

ZNF141 siRNA (h): sc-89148

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which is thought to interact with KAP1, thereby recruiting histone modifying proteins. As a member of the Krüppel C₂H₂-type zinc-finger protein family, ZNF141 (zinc finger protein 141) is a 474 amino acid nuclear protein that contains one KRAB domain and 11 C₂H₂-type zinc fingers. ZNF141 acts as a transcription repressor and is expressed ubiquitously at low levels. The gene encoding ZNF141 maps within a region of human chromosome 4p16.3 that is deleted in Wolf-Hirschhorn syndrome, a condition that leads to severe abnormalities such as microencephaly, mental retardation, cleft lip and/or cleft palate, seizures and poor muscle tone.

REFERENCES

1. Wolf, U., Reinwein, H., Porsch, R., Schröter, R. and Baitsch, H. 1965. Deficiency on the short arms of a chromosome No. 4. *Humangenetik* 1: 397-413.
2. Hirschhorn, K., Cooper, H.L. and Firschein, I.L. 1965. Deletion of short arms of chromosome 4-5 in a child with defects of midline fusion. *Humangenetik* 1: 479-482.
3. Bakker, H.R. and Waring, H. 1976. Experimental diabetes insipidus in a marsupial, *Macropus eugenii* (Desmarest). *J. Endocrinol.* 69: 149-157.
4. Overbeke, A.J. 1989. Donor blood or autologous blood, again the question. *Ned. Tijdschr. Geneesk.* 133: 1156-1157.
5. Tommerup, N., Aagaard, L., Lund, C.L., Boel, E., Baxendale, S., Bates, G.P., Lehrach, H. and Vissing, H. 1993. A zinc-finger gene ZNF141 mapping at 4p16.3/D4S90 is a candidate gene for the Wolf-Hirschhorn (4p-) syndrome. *Hum. Mol. Genet.* 2: 1571-1575.
6. Krishna, S.S., Majumdar, I. and Grishin, N.V. 2003. Structural classification of zinc fingers: survey and summary. *Nucleic Acids Res.* 31: 532-550.
7. Mihara, M., Yoshida, Y., Tsukamoto, T., Inada, K., Nakanishi, Y., Yagi, Y., Imai, K., Sugimura, T., Tatematsu, M. and Ushijima, T. 2006. Methylation of multiple genes in gastric glands with intestinal metaplasia: a disorder with polyclonal origins. *Am. J. Pathol.* 169: 1643-1651.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2010. Johns Hopkins University, Baltimore, MD. MIM Number: 194648. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: ZNF141 (human) mapping to 4p16.3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

ZNF141 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ZNF141 shRNA Plasmid (h): sc-89148-SH and ZNF141 shRNA (h) Lentiviral Particles: sc-89148-V as alternate gene silencing products.

For independent verification of ZNF141 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-89148A, sc-89148B and sc-89148C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ZNF141 siRNA (h) is recommended for the inhibition of ZNF141 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ZNF141 gene expression knockdown using RT-PCR Primer: ZNF141 (h)-PR: sc-89148-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.