

# UGT2B15 siRNA (h): sc-89282

## BACKGROUND

UDP-glucuronosyltransferase isoenzymes (UGTs) catalyze the glucuronidation of small lipophilic molecules, which regulates the bioactivity and metabolic fate of a wide range of endo- and xenobiotics. Glucuronidation increases the polarity of lipophilic molecules and facilitates their entry into aqueous compartments and their ultimate excretion. In essence, glucuronidation provides a protective function by terminating or attenuating the biological activity of its substrates. The UGT2B family of isoenzymes are highly expressed in liver, but are also detected in several non-hepatic tissues, including skin, breast, prostate, intestine, placenta, and lung. Therefore, the UGT2B family may preferentially modulate steroid metabolism and excretion in addition to bile acids and xenobiotics. The human UGT2B genes localize as a cluster on chromosome 4q13.2. UGT2B15 (UDP-glucuronosyltransferase 2B15), also known as HLUG4 and UDP-glucuronosyltransferase 2B8, is a 530 amino acid single-pass membrane protein that is a member of the UGT2B family.

## REFERENCES

1. Coffman, B.L., et al. 1990. Characterization and primary sequence of a human hepatic microsomal estradiol UDP glucuronosyltransferase. *Arch. Biochem. Biophys.* 281: 170-175.
2. Chen, F., et al. 1993. Characterization of a cloned human dihydrotestosterone/androstenediol UDP-glucuronosyltransferase and its comparison to other steroid isoforms. *Biochemistry* 32: 10648-10657.
3. Green, M.D., et al. 1994. Stable expression of a human liver UDP-glucuronosyltransferase (UGT2B15) with activity toward steroid and xenobiotic substrates. *Drug Metab. Dispos.* 22: 799-805.
4. Monaghan, G., et al. 1994. Isolation of a human YAC contig encompassing a cluster of UGT2 genes and its regional localization to chromosome 4q13. *Genomics* 23: 496-499.
5. Levesque, E., et al. 1997. Isolation and characterization of UGT2B15(Y85): a UDP-glucuronosyltransferase encoded by a polymorphic gene. *Pharmacogenetics* 7: 317-325.
6. Toide, K., et al. 2002. A major genotype in UDP-glucuronosyltransferase 2B15. *Drug Metab. Pharmacokinet.* 17: 164-166.
7. Iida, A., et al. 2002. Catalog of 86 single-nucleotide polymorphisms (SNPs) in three uridine diphosphate glycosyltransferase genes: UGT2A1, UGT2B15, and UGT8. *J. Hum. Genet.* 47: 505-510.
8. Mackenzie, P.I., et al. 2005. Nomenclature update for the mammalian UDP glycosyltransferase (UGT) gene superfamily. *Pharmacogenet. Genomics* 15: 677-685.
9. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2009. Johns Hopkins University, Baltimore, MD. MIM Number: 600069. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: UGT2B15 (human) mapping to 4q13.2.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

UGT2B15 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see UGT2B15 shRNA Plasmid (h): sc-89282-SH and UGT2B15 shRNA (h) Lentiviral Particles: sc-89282-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

UGT2B15 siRNA (h) is recommended for the inhibition of UGT2B15 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor UGT2B15 gene expression knockdown using RT-PCR Primer: UGT2B15 (h)-PR: sc-89282-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.