



Polycystin-1L1 siRNA (h): sc-89442

BACKGROUND

Polycystin-1L1, also known as PKD1L1 (polycystic kidney disease 1 like 1), PRO19563 or PC1-like 1 protein, is a 2,849 amino acid multi-pass membrane protein. Belonging to the Polycystin family, Polycystin-1L1 contains 11 transmembrane domains, one GPS domain, two Ig-like PKD domains, one LH2/polycystin-1, lipoxigenase, alpha-toxin (PLAT) domain, one small receptor for egg jelly (REJ) domain and one coiled-coil domain. Encoded by a gene that maps to human chromosome 7p12.3, Polycystin-1L1 contains 58 exons and exists as two alternatively spliced isoforms. Polycystin-1L1 shares significant homology with all known Polycystins, and 61% sequence identity with its mouse homolog. Polycystin-1L1 is expressed in fetal and adult heart and in testis. Polycystin-1L1 is also expressed strongly in Leydig cells, a testosterone production source, and may play a role in the male reproductive system. Defects in the gene encoding Polycystin-1L1 may be associated with Polycystic kidney disease, a progressive disorder characterized by the presence of cysts in the kidneys.

REFERENCES

1. Yuasa, T., Venugopal, B., Weremowicz, S., Morton, C.C., Guo, L. and Zhou, J. 2002. The sequence, expression, and chromosomal localization of a novel polycystic kidney disease 1-like gene, PKD1L1, in human. *Genomics* 79: 376-386.
2. Luo, Y., Vassilev, P.M., Li, X., Kawanabe, Y. and Zhou, J. 2003. Native polycystin 2 functions as a plasma membrane Ca^{2+} -permeable cation channel in renal epithelia. *Mol. Cell. Biol.* 23: 2600-2607.
3. Clapham, D.E. 2003. TRP channels as cellular sensors. *Nature* 426: 517-524.
4. Lakkis, M. and Zhou, J. 2003. Molecular complexes formed with polycystins. *Nephron Exp. Nephrol.* 93: e3-e8.
5. Nauli, S.M. and Zhou, J. 2004. Polycystins and mechanosensation in renal and nodal cilia. *Bioessays* 26: 844-856.
6. Yuasa, T., Takakura, A., Denker, B.M., Venugopal, B. and Zhou, J. 2004. Polycystin-1L2 is a novel G protein-binding protein. *Genomics* 84: 126-138.
7. Zhou, J. 2009. Polycystins and primary cilia: primers for cell cycle progression. *Annu. Rev. Physiol.* 71: 83-113.

CHROMOSOMAL LOCATION

Genetic locus: PKD1L1 (human) mapping to 7p12.3.

PRODUCT

Polycystin-1L1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Polycystin-1L1 shRNA Plasmid (h): sc-89442-SH and Polycystin-1L1 shRNA (h) Lentiviral Particles: sc-89442-V as alternate gene silencing products.

For independent verification of Polycystin-1L1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-89442A, sc-89442B and sc-89442C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Polycystin-1L1 siRNA (h) is recommended for the inhibition of Polycystin-1L1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Polycystin-1L1 gene expression knockdown using RT-PCR Primer: Polycystin-1L1 (h)-PR: sc-89442-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.