

ATP5J2 siRNA (h): sc-89543

BACKGROUND

Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. It is composed of two linked multi-subunit complexes: the soluble catalytic core, F₁, and the membrane-spanning component, F₀, which comprises the proton channel. ATP5J2, also known as ATP synthase subunit f, mitochondrial, is a 94 amino acid mitochondrion inner membrane that belongs to the ATPase F chain family. Mitochondrial dysfunction is prominent in Alzheimer's disease (AD). A failure of one or more of the mitochondrial electron transport chain enzymes, or of F₁F₀-ATPase (ATP synthase), could compromise brain energy stores, generate damaging reactive oxygen species (ROS), and lead to neuronal death. Existing as two alternatively spliced isoforms, the ATP5J2 gene is conserved in chimpanzee, dog, canine, bovine, rat, chicken and zebrafish, and maps to human chromosome 7q22.1.

REFERENCES

1. Elston, T., Wang, H. and Oster, G. 1998. Energy transduction in ATP synthase. *Nature* 391: 510-513.
2. Wang, H. and Oster, G. 1998. Energy transduction in the F₁ motor of ATP synthase. *Nature* 396: 279-282.
3. Kinoshita, K., Yasuda, R. and Noji, H. 2000. F₁-ATPase: a highly efficient rotary ATP machine. *Essays Biochem.* 35: 3-18.
4. Aggeler, R., Coons, J., Taylor, S.W., Ghosh, S.S., Garcia, J.J., Capaldi, R.A. and Marusich, M.F. 2002. A functionally active human F₁F₀ ATPase can be purified by immunocapture from heart tissue and fibroblast cell lines. Subunit structure and activity studies. *J. Biol. Chem.* 277: 33906-33912.
5. Bosetti, F., Brizzi, F., Barogi, S., Mancuso, M., Siciliano, G., Tendi, E.A., Murri, L., Rapoport, S.I. and Solaini, G. 2002. Cytochrome c oxidase and mitochondrial F₁F₀-ATPase (ATP synthase) activities in platelets and brain from patients with Alzheimer's disease. *Neurobiol. Aging* 23: 371-376.
6. Leyva, J.A., Bianchet, M.A. and Amzel, L.M. 2003. Understanding ATP synthesis: structure and mechanism of the F₁-ATPase (review). *Mol. Membr. Biol.* 20: 27-33.
7. Scherer, S.W., Cheung, J., MacDonald, J.R., Osborne, L.R., Nakabayashi, K., Herbrick, J.A., Carson, A.R., Parker-Katiraei, L., Skaug, J., Khaja, R., Zhang, J., Hudek, A.K., Li, M., Haddad, M., Duggan, G.E., et al. 2003. Human chromosome 7: DNA sequence and biology. *Science* 300: 767-772.
8. Cross, R.L. 2004. Molecular motors: turning the ATP motor. *Nature* 427: 407-408.

CHROMOSOMAL LOCATION

Genetic locus: ATP5J2 (human) mapping to 7q22.1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

ATP5J2 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ATP5J2 shRNA Plasmid (h): sc-89543-SH and ATP5J2 shRNA (h) Lentiviral Particles: sc-89543-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ATP5J2 siRNA (h) is recommended for the inhibition of ATP5J2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.