



# hyccin siRNA (h): sc-89596

## BACKGROUND

Hyccin, also known as FAM126A, HCC, HLD5, HYCC1 or DRCTNNB1A, is a 521 amino acid cytoplasmic protein that is widely expressed with highest levels found in heart, brain, placenta, spleen and testis. Belonging to the FAM126 family, hyccin may play a role in the  $\beta$ -catenin/Lef signaling pathway. Hyccin is likely involved in the process of myelination of the central and peripheral nervous system. Defects in the gene encoding hyccin are the cause of leukodystrophy hypomyelinating type 5 (HLD5), which is characterized by congenital cataract, progressive neurologic impairment and diffuse myelin deficiency. Individuals affected by HLD5 experience progressive pyramidal and cerebellar dysfunction along with muscle weakness in the lower limbs. Hyccin exists as two alternatively spliced isoforms and is encoded by a gene located on human chromosome 7.

## REFERENCES

1. Kawasoe, T., Furukawa, Y., Daigo, Y., Nishiwaki, T., Ishiguro, H., Fujita, M., Satoh, S., Miwa, N., Nagasawa, Y., Miyoshi, Y., Ogawa, M. and Nakamura, Y. 2000. Isolation and characterization of a novel human gene, DRCTNNB1A, the expression of which is down-regulated by  $\beta$ -catenin. *Cancer Res.* 60: 3354-3358.
2. Zara, F., Biancheri, R., Bruno, C., Bordo, L., Assereto, S., Gazzero, E., Sotgia, F., Wang, X.B., Gianotti, S., Stringara, S., Pedemonte, M., Uziel, G., Rossi, A., Schenone, A., Tortori-Donati, P., van der Knaap, M.S., et al. 2006. Deficiency of hyccin, a newly identified membrane protein, causes hypomyelination and congenital cataract. *Nat. Genet.* 38: 1111-1113.
3. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610531. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Biancheri, R., Zara, F., Bruno, C., Rossi, A., Bordo, L., Gazzero, E., Sotgia, F., Pedemonte, M., Scapolan, S., Bado, M., Uziel, G., Bugiani, M., Lamba, L.D., Costa, V., Schenone, A., Rozemuller, A.J., Tortori-Donati, P., et al. 2007. Phenotypic characterization of hypomyelination and congenital cataract. *Ann. Neurol.* 62: 121-127.
5. Rossi, A., Biancheri, R., Zara, F., Bruno, C., Uziel, G., van der Knaap, M.S., Minetti, C. and Tortori-Donati, P. 2008. Hypomyelination and congenital cataract: neuroimaging features of a novel inherited white matter disorder. *AJNR Am. J. Neuroradiol.* 29: 301-305.
6. Bizzi, A., Castelli, G., Bugiani, M., Barker, P.B., Herskovits, E.H., Danesi, U., Erbetta, A., Moroni, I., Farina, L. and Uziel, G. 2008. Classification of childhood white matter disorders using proton MR spectroscopic imaging. *AJNR Am. J. Neuroradiol.* 29: 1270-1275.
7. Ugur, S.A. and Tolun, A. 2008. A deletion in DRCTNNB1A associated with hypomyelination and juvenile onset cataract. *Eur. J. Hum. Genet.* 16: 261-264.

## CHROMOSOMAL LOCATION

Genetic locus: FAM126A (human) mapping to 7p15.3.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

hyccin siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see hyccin shRNA Plasmid (h): sc-89596-SH and hyccin shRNA (h) Lentiviral Particles: sc-89596-V as alternate gene silencing products.

For independent verification of hyccin (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-89596A, sc-89596B and sc-89596C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

hyccin siRNA (h) is recommended for the inhibition of hyccin expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor hyccin gene expression knockdown using RT-PCR Primer: hyccin (h)-PR: sc-89596-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.