

ATP10A siRNA (h): sc-89984

BACKGROUND

ATPase class V type 10A (ATP10A), also known as aminophospholipid translocase VA (ATPVA) or probable phospholipid-transporting ATPase VA, is a 1,499 amino acid protein belonging to the cation transport ATPase (P-type) family. Localized to the cell membrane, ATP10A is widely expressed in various tissues, with highest levels in kidney, lung, brain, prostate, testis, ovary and small intestine. ATP10A transports phosphatidylserine and phosphatidylethanolamine from one side of a membrane lipid bilayer to another. The gene encoding ATP10A is an imprinted gene that is maternally expressed. Defects in this gene lead to Angelman syndrome (AS), also known as happy puppet syndrome. AS is characterized by mental retardation, movement or balance disorder, characteristic abnormal behaviors and severe limitations in speech and language.

REFERENCES

1. Kayashima, T., et al. 2003. Atp10a, the mouse ortholog of the human imprinted ATP10A gene, escapes genomic imprinting. *Genomics* 81: 644-647.
2. Kayashima, T., et al. 2003. On the conflicting reports of imprinting status of mouse ATP10a in the adult brain: strain-background-dependent imprinting? *J. Hum. Genet.* 48: 492-493
3. Online Mendelian Inheritance in Man, OMIM[™]. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 605855. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Kim, K.P., et al. 2007. Gene-specific vulnerability to imprinting variability in human embryonic stem cell lines. *Genome Res.* 17: 1731-1742.
5. Kato, C., et al. 2008. Association study of the 15q11-q13 maternal expression domain in Japanese autistic patients. *Am. J. Med. Genet. B Neuropsychiatr. Genet.* 147B: 1008-1012.
6. Giardina, E., et al. 2008. A multiplex molecular assay for the detection of uniparental disomy for human chromosome 15q11.2. *Electrophoresis* 29: 4775-4779.

CHROMOSOMAL LOCATION

Genetic locus: ATP10A (human) mapping to 15q12.

PRODUCT

ATP10A siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ATP10A shRNA Plasmid (h): sc-89984-SH and ATP10A shRNA (h) Lentiviral Particles: sc-89984-V as alternate gene silencing products.

For independent verification of ATP10A (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-89984A, sc-89984B and sc-89984C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ATP10A siRNA (h) is recommended for the inhibition of ATP10A expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

ATP10A (G-9): sc-514650 is recommended as a control antibody for monitoring of ATP10A gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ATP10A gene expression knockdown using RT-PCR Primer: ATP10A (h)-PR: sc-89984-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.