

IDI1 siRNA (h): sc-90838

BACKGROUND

IDI1 (isopentenyl-diphosphate δ isomerase 1), also known as IPP1 or IPP11, is a 227 amino acid member of the IPP isomerase type I family and is involved in cholesterol biosynthesis. Localized to the peroxisome, IDI1 catalytically converts isopentenyl diphosphate (IPP) to its electrophilic isomer, dimethylallyl diphosphate (DMAPP). Specifically, IDI1 uses magnesium as a cofactor to catalyze the 1,3-allylic rearrangement of IPP, thus creating DMAPP, a substrate for subsequent reactions that synthesize farnesyl diphosphate and, ultimately, cholesterol. Defects in the gene encoding IDI1 may be associated with peroxisomal deficiency diseases, such as Zellweger syndrome, a congenital disorder caused by a reduction in the number of peroxisomes. Individuals affected with this disorder generally exhibit lack of muscle tone, an enlarged liver, mental retardation and, in some cases, death.

REFERENCES

1. Xuan, J.W., et al. 1994. A human promyelocyte mRNA transiently induced by TPA is homologous to yeast IPP isomerase. *Genomics* 20: 129-131.
2. Krisans, S.K., et al. 1994. Farnesyl-diphosphate synthase is localized in peroxisomes. *J. Biol. Chem.* 269: 14165-14169.
3. Hahn, F.M., et al. 1996. Human isopentenyl diphosphate: dimethylallyl diphosphate isomerase: overproduction, purification, and characterization. *Arch. Biochem. Biophys.* 332: 30-34.
4. Paton, V.G., et al. 1997. Cloning and subcellular localization of hamster and rat isopentenyl diphosphate dimethylallyl diphosphate isomerase. A PTS1 motif targets the enzyme to peroxisomes. *J. Biol. Chem.* 272: 18945-18950.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604055. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Breitling, R., et al. 2003. Isopentenyl-diphosphate isomerases in human and mouse: evolutionary analysis of a mammalian gene duplication. *J. Mol. Evol.* 57: 282-291.
7. Clizbe, D.B., et al. 2007. IDI2, a second isopentenyl diphosphate isomerase in mammals. *J. Biol. Chem.* 282: 6668-6676.

CHROMOSOMAL LOCATION

Genetic locus: IDI1 (human) mapping to 10p15.3.

PRODUCT

IDI1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see IDI1 shRNA Plasmid (h): sc-90838-SH and IDI1 shRNA (h) Lentiviral Particles: sc-90838-V as alternate gene silencing products.

For independent verification of IDI1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-90838A, sc-90838B and sc-90838C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

IDI1 siRNA (h) is recommended for the inhibition of IDI1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

IDI1 (XY-7): sc-100550 is recommended as a control antibody for monitoring of IDI1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor IDI1 gene expression knockdown using RT-PCR Primer: IDI1 (h)-PR: sc-90838-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.