

# NAP1L3 siRNA (h): sc-91058

## BACKGROUND

Proper nucleosome assembly is critical for compacting DNA into chromatin. NAP1 (nucleosome assembly protein 1) is a nuclear protein that acts as a transcriptional regulator and functions in nucleosome assembly. NAP1L3 (nucleosome assembly protein 1-like 3), also known as MB20 or NPL3, is a 506 amino acid nuclear protein belonging to the nucleosome assembly protein (NAP) family. Expressed in human brain with weak expression in heart, NAP1L3 is encoded by a gene mapping to human chromosome Xq21.32, which is in close proximity to a region closely linked to several X-linked mental retardation syndromes. Containing nearly 153 million base pairs and housing over 1,000 genes, chromosome X acts in conjunction with chromosome Y to determine sex. There are a number of conditions related to an abnormal number and combination of sex chromosomes, some of which include Turner's syndrome, color blindness, hemophilia and Duchenne muscular dystrophy.

## REFERENCES

1. Watanabe, T.K., Fujiwara, T., Nakamura, Y., Hirai, Y., Maekawa, H. and Takahashi, E. 1996. Cloning, expression pattern and mapping to Xq of NAP1L3, a gene encoding a peptide homologous to human and yeast nucleosome assembly proteins. *Cytogenet. Cell Genet.* 74: 281-285.
2. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 300117. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Uprichard, J. and Perry, D.J. 2002. Factor X deficiency. *Blood Rev.* 16: 97-110.
4. Sharpe, L.T., de Luca, E., Hansen, T., Jägle, H. and Gegenfurtner, K.R. 2006. Advantages and disadvantages of human dichromacy. *J. Vis.* 6: 213-223.
5. Helderma-van den Enden, A.T., de Jong, R., den Dunnen, J.T., Houwing-Duistermaat, J.J., Kneppers, A.L., Ginjaar, H.B., Breuning, M.H. and Bakker, E. 2009. Recurrence risk due to germ line mosaicism: Duchenne and Becker muscular dystrophy. *Clin. Genet.* 75: 465-472.
6. Mullaney, R. and Murphy, D. 2009. Turner syndrome: neuroimaging findings: structural and functional. *Dev. Disabil. Res. Rev.* 15: 279-283.
7. Makishima, T., King, K., Brewer, C.C., Zalewski, C.K., Butman, J., Bakalov, V.K., Bondy, C. and Griffith, A.J. 2009. Otolaryngologic markers for the early diagnosis of Turner syndrome. *Int. J. Pediatr. Otorhinolaryngol.* 73: 1564-1567.

## CHROMOSOMAL LOCATION

Genetic locus: NAP1L3 (human) mapping to Xq21.32.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

NAP1L3 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see NAP1L3 shRNA Plasmid (h): sc-91058-SH and NAP1L3 shRNA (h) Lentiviral Particles: sc-91058-V as alternate gene silencing products.

For independent verification of NAP1L3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-91058A, sc-91058B and sc-91058C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

NAP1L3 siRNA (h) is recommended for the inhibition of NAP1L3 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor NAP1L3 gene expression knockdown using RT-PCR Primer: NAP1L3 (h)-PR: sc-91058-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.