

MAGE-H1 siRNA (h): sc-91276

BACKGROUND

The melanoma-associated antigen (MAGE) family consists of a number of antigens recognized by cytotoxic T lymphocytes. The MAGE genes were initially isolated from different kinds of tumors and, based on their virtually exclusive tumor-specific expression in adult tissues, they have been used as targets for cancer immunotherapy. MAGE genes encode for tumor-rejection antigens that are expressed in tumors of different histologic types and in normal testis and placenta. MAGE-H1 (Melanoma-associated antigen H1), also known as Restin or APR1 (Apoptosis-related protein 1), is a 219 amino acid protein that contains a type II MAGE homology domain (MHD). Enhanced ligand stimulation promotes MAGE-H1 interaction with the type II death domain of NGFR p75. It is suggested that MAGE-H1 accelerates differentiation in response to nerve growth factor in cells.

REFERENCES

1. De Plaen, E., Arden, K., Traversari, C., Gaforio, J.J., Szikora, J.P., De Smet, C., Brasseur, F., van der Bruggen, P., Lethé, B. and Lurquin, C. 1994. Structure, chromosomal localization, and expression of 12 genes of the MAGE family. *Immunogenetics* 40: 360-369.
2. Lucas, S., Brasseur, F. and Boon, T. 1999. A new MAGE gene with ubiquitous expression does not code for known MAGE antigens recognized by T cells. *Cancer Res.* 59: 4100-4103.
3. Serrano, A., Lethé, B., Delroisse, J.M., Lurquin, C., De Plaen, E., Brasseur, F., Rimoldi, D. and Boon, T. 1999. Quantitative evaluation of the expression of MAGE genes in tumors by limiting dilution of cDNA libraries. *Int. J. Cancer* 83: 664-669.
4. Kobayashi, Y., Higashi, T., Nouse, K., Nakatsukasa, H., Ishizaki, M., Kaneyoshi, T., Toshikuni, N., Kariyama, K., Nakayama, E. and Tsuji, T. 2000. Expression of MAGE, GAGE and BAGE genes in human liver diseases: utility as molecular markers for hepatocellular carcinoma. *J. Hepatol.* 32: 612-617.
5. Tcherpakov, M., Bronfman, F.C., Conticello, S.G., Vaskovsky, A., Levy, Z., Niinobe, M., Yoshikawa, K., Arenas, E. and Fainzilber, M. 2002. The p75 neurotrophin receptor interacts with multiple MAGE proteins. *J. Biol. Chem.* 277: 49101-49104.
6. Barker, P.A. and Salehi, A. 2002. The MAGE proteins: emerging roles in cell cycle progression, apoptosis, and neurogenetic disease. *J. Neurosci. Res.* 67: 705-712.
7. Shao, J.B. and Chen, Z. 2003. Expression of MAGE, GAGE, and BAGE genes in human hepatocellular carcinoma. *Zhonghua Gan Zang Bing Za Zhi* 11: 142-144.
8. Albrecht, D.E. and Froehner, S.C. 2004. DAMAGE, a novel α -dystrobrevin-associated MAGE protein in dystrophin complexes. *J. Biol. Chem.* 279: 7014-7023.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: MAGEH1 (human) mapping to Xp11.21.

PRODUCT

MAGE-H1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see MAGE-H1 shRNA Plasmid (h): sc-91276-SH and MAGE-H1 shRNA (h) Lentiviral Particles: sc-91276-V as alternate gene silencing products.

For independent verification of MAGE-H1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-91276A, sc-91276B and sc-91276C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

MAGE-H1 siRNA (h) is recommended for the inhibition of MAGE-H1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor MAGE-H1 gene expression knockdown using RT-PCR Primer: MAGE-H1 (h)-PR: sc-91276-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.