

GASP-2 siRNA (h): sc-91285

BACKGROUND

G protein-coupled receptors play a role in many different stimulus-response pathways. G protein-coupled receptors mediate extracellular signals into intracellular signals (G protein activation). They respond to a wide variety of signaling molecules, including hormones, neurotransmitters and other proteins and peptides. GASP-2 (G protein-coupled receptor associated sorting protein 2), also known as GPRASP2, is an 838 amino acid protein that regulates a number of G protein-coupled receptors, such as CT-R (calcitonin receptor) and mAChR M1 (muscarinic acetylcholine receptor M1), through interactions with their cytoplasmic tails. Expressed primarily in brain, GASP-2 is a member of the GPRASP family and forms a complex with Huntingtin, with which it is thought to influence receptor trafficking.

REFERENCES

1. Lee, D.K., George, S.R. and O'Dowd, B.F. 2002. Novel G protein-coupled receptor genes expressed in the brain: continued discovery of important therapeutic targets. *Expert Opin. Ther. Targets* 6: 185-202
2. Simonin, F., Karcher, P., Boeuf, J.J., Matifas, A. and Kieffer, B.L. 2004. Identification of a novel family of G protein-coupled receptor associated sorting proteins. *J. Neurochem.* 89: 766-775.
3. Horn, S.C., Lalowski, M., Goehler, H., Dröge, A., Wanker, E.E. and Stelzl, U. 2006. Huntingtin interacts with the receptor sorting family protein GASP-2. *J. Neural. Transm.* 113: 1081-1090.
4. Rozenfeld, R. and Devi, L.A. 2010. Exploring a role for heteromerization in GPCR signalling specificity. *Biochem. J.* 433: 11-18.
5. Costanzi, S. 2010. Modeling G protein-coupled receptors: a concrete possibility. *Chim. Oggi* 28: 26-31.
6. Borroto-Escuela, D.O., Correia, P.A., Romero-Fernandez, W., Narvaez, M., Fuxe, K., Ciruela, F. and Garriga, P. 2011. Muscarinic receptor family interacting proteins: role in receptor function. *J. Neurosci. Methods* 195: 161-169.

CHROMOSOMAL LOCATION

Genetic locus: GPRASP2 (human) mapping to Xq22.1.

PRODUCT

GASP-2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GASP-2 shRNA Plasmid (h): sc-91285-SH and GASP-2 shRNA (h) Lentiviral Particles: sc-91285-V as alternate gene silencing products.

For independent verification of GASP-2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-91285A, sc-91285B and sc-91285C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GASP-2 siRNA (h) is recommended for the inhibition of GASP-2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

GASP-2 (F-6): sc-390664 is recommended as a control antibody for monitoring of GASP-2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GASP-2 gene expression knockdown using RT-PCR Primer: GASP-2 (h)-PR: sc-91285-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.