

P2RY8 siRNA (h): sc-91574

BACKGROUND

G protein-coupled receptors (GPRs), also known as seven transmembrane receptors, heptahelical receptors or 7TM receptors, comprise a superfamily of proteins that play a role in many different stimulus-response pathways. G protein-coupled receptors translate extracellular signals into intracellular signals (G protein activation) and they respond to a variety of signaling molecules, such as hormones and neurotransmitters. P2RY8 (purinergic receptor P2Y, G protein-coupled, 8), also known as P2Y8, is a 359 amino acid multi-pass membrane protein that localizes to the cell membrane and belongs to the G protein-coupled receptor family. Expressed at low levels in lung, heart and kidney, P2RY8 may function as a receptor for purines that are coupled to G proteins and may also play a role in mental retardation.

REFERENCES

1. Adrian, K., Bernhard, M.K., Breiting, H.G. and Ogilvie, A. 2000. Expression of purinergic receptors (ionotropic P2X1-7 and metabotropic P2Y1-11) during myeloid differentiation of HL-60 cells. *Biochim. Biophys. Acta* 1492: 127-138.
2. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 300525. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Cantagrel, V., Lossi, A.M., Boulanger, S., Depetris, D., Mattei, M.G., Gecz, J., Schwartz, C.E., Van Maldergem, L. and Villard, L. 2004. Disruption of a new X linked gene highly expressed in brain in a family with two mentally retarded males. *J. Med. Genet.* 41: 736-742.
4. Fujiwara, S., Yamashita, Y., Choi, Y.L., Watanabe, H., Kurashina, K., Soda, M., Enomoto, M., Hatanaka, H., Takada, S., Ozawa, K. and Mano, H. 2007. Transforming activity of purinergic receptor P2Y, G protein-coupled, 8 revealed by retroviral expression screening. *Leuk. Lymphoma* 48: 978-986.
5. Storlazzi, C.T., Albano, F., Lo Cunsolo, C., Doglioni, C., Guastadisegni, M.C., Impera, L., Lonoce, A., Funes, S., Macrì, E., Iuzzolino, P., Panagopoulos, I., Specchia, G. and Rocchi, M. 2007. Upregulation of the Sox-5 by promoter swapping with the P2RY8 gene in primary splenic follicular lymphoma. *Leukemia* 21: 2221-2225.
6. Regard, J.B., Sato, I.T. and Coughlin, S.R. 2008. Anatomical profiling of G protein-coupled receptor expression. *Cell* 135: 561-571.

CHROMOSOMAL LOCATION

Genetic locus: P2RY8 (human) mapping to Xp22.33/Yp11.32.

PRODUCT

P2RY8 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see P2RY8 shRNA Plasmid (h): sc-91574-SH and P2RY8 shRNA (h) Lentiviral Particles: sc-91574-V as alternate gene silencing products.

For independent verification of P2RY8 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-91574A, sc-91574B and sc-91574C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCL, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

P2RY8 siRNA (h) is recommended for the inhibition of P2RY8 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor P2RY8 gene expression knockdown using RT-PCR Primer: P2RY8 (h)-PR: sc-91574-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.