



# ARL15 siRNA (h): sc-91713

## BACKGROUND

ADP-ribosylation factors (ARFs) are highly conserved guanine nucleotide binding proteins that enhance the ADP-ribosyltransferase activity of cholera toxin. ARFs are important in eukaryotic vesicular trafficking pathways and play an essential role in the activation of phospholipase D (PC-PLD). ARL15 (ADP-ribosylation factor-like 15), also known as ARFRP2 (ARF-related protein 2), is a 204 amino acid protein that belongs to the small GTPase superfamily and the ARF family. Conserved in chimpanzee, canine, bovine, mouse, rat, chicken and zebrafish, ARL15 participates in both GTP and nucleotide binding. ARL15 variants influence levels of Acip30, an adipocyte-derived protein that is particularly heritable and inversely associated with risk of type 2 diabetes mellitus and coronary heart disease. ARL15 is encoded by a gene that maps to human chromosome 5q11.2.

## REFERENCES

1. Shimada, T., et al. 1983. Induction of transformation by six classes of chemical carcinogens in adult rat liver epithelial cells. *Cancer Res.* 43: 5087-5092.
2. Furukawa, K., et al. 1987. Enrichment and characterization of clonogenic epithelial cells from adult rat liver and initiation of epithelial cell strains. *In Vitro Cell. Dev. Biol.* 23: 339-348.
3. Somani, S.M., et al. 1990. Eseroline, a metabolite of physostigmine, induces neuronal cell death. *Toxicol. Appl. Pharmacol.* 106: 28-37.
4. Gillingham, A.K., et al. 2007. The small G proteins of the Arf family and their regulators. *Annu. Rev. Cell Dev. Biol.* 23: 579-611.
5. Macqueen, D.J., et al. 2008. Evolution of follistatin in teleosts revealed through phylogenetic, genomic and expression analyses. *Dev. Genes Evol.* 218: 1-14.
6. Richards, J.B., et al. 2009. A genome-wide association study reveals variants in ARL15 that influence adiponectin levels. *PLoS Genet.* 5: e1000768.
7. Barton, V.N., et al. 2010. Unique molecular characteristics of pediatric myxopapillary ependymoma. *Brain Pathol.* 20: 560-570.
8. Cook, J.R., et al. 2010. Hypoadiponectinemia—cause or consequence of human “insulin resistance?” *J. Clin. Endocrinol. Metab.* 95: 1544-1554.
9. Glessner, J.T., et al. 2010. A genome-wide study reveals copy number variants exclusive to childhood obesity cases. *Am. J. Hum. Genet.* 87: 661-666.

## CHROMOSOMAL LOCATION

Genetic locus: ARL15 (human) mapping to 5q11.2.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

ARL15 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ARL15 shRNA Plasmid (h): sc-91713-SH and ARL15 shRNA (h) Lentiviral Particles: sc-91713-V as alternate gene silencing products.

For independent verification of ARL15 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-91713A, sc-91713B and sc-91713C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ARL15 siRNA (h) is recommended for the inhibition of ARL15 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ARL15 gene expression knockdown using RT-PCR Primer: ARL15 (h)-PR: sc-91713-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.