

# BSPRY siRNA (h): sc-92522

## BACKGROUND

BSPRY (B-box and SPRY domain containing) is a 402 amino acid protein that localizes to both the membrane and the cytoplasm and contains one B box-type zinc finger and one B30.2/SPRY domain. Existing as two alternatively spliced isoforms, BSPRY interacts with TRPV5 and TRPV6 and is thought to regulate the transport of calcium across the epithelium, probably by inhibiting the activity of TRPV proteins. The gene encoding BSPRY maps to human chromosome 9, which houses over 900 genes and comprises nearly 4% of the human genome. Hereditary hemorrhagic telangiectasia, which is characterized by harmful vascular defects, and Familial dysautonomia, are both associated with chromosome 9. Notably, chromosome 9 encompasses the largest interferon family gene cluster.

## REFERENCES

1. Schenker, T. and Trueb, B. 2000. BSPRY, a novel protein of the Ro-Ret family. *Biochim. Biophys. Acta* 1493: 255-258.
2. Gisler, S.M., Stagljar, I., Traebert, M., Bacic, D., Biber, J. and Murer, H. 2001. Interaction of the type IIa Na/P<sub>i</sub> co-transporter with PDZ proteins. *J. Biol. Chem.* 276: 9206-9213.
3. Birkenfeld, J., Kartmann, B., Anliker, B., Ono, K., Schlötcke, B., Betz, H. and Roth, D. 2003. Characterization of zetin 1/rBSPRY, a novel binding partner of 14-3-3 proteins. *Biochem. Biophys. Res. Commun.* 302: 526-533.
4. van Abel, M., Hoenderop, J.G. and Bindels, R.J. 2005. The epithelial calcium channels TRPV5 and TRPV6: regulation and implications for disease. *Naunyn Schmiedeberg's Arch. Pharmacol.* 371: 295-306.
5. Nijenhuis, T., Hoenderop, J.G. and Bindels, R.J. 2005. TRPV5 and TRPV6 in Ca<sup>2+</sup> (re)absorption: regulating Ca<sup>2+</sup> entry at the gate. *Pflugers Arch.* 451: 181-192.
6. van de Graaf, S.F., van der Kemp, A.W., van den Berg, D., van Oorschot, M., Hoenderop, J.G. and Bindels, R.J. 2006. Identification of BSPRY as a novel auxiliary protein inhibiting TRPV5 activity. *J. Am. Soc. Nephrol.* 17: 26-30.
7. Gkika, D., Hsu, Y.J., van der Kemp, A.W., Christakos, S., Bindels, R.J. and Hoenderop, J.G. 2006. Critical role of the epithelial Ca<sup>2+</sup> channel TRPV5 in active Ca<sup>2+</sup> reabsorption as revealed by TRPV5/calbindin-D28K knockout mice. *J. Am. Soc. Nephrol.* 17: 3020-3027.

## CHROMOSOMAL LOCATION

Genetic locus: BSPRY (human) mapping to 9q32.

## PRODUCT

BSPRY siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see BSPRY shRNA Plasmid (h): sc-92522-SH and BSPRY shRNA (h) Lentiviral Particles: sc-92522-V as alternate gene silencing products.

For independent verification of BSPRY (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-92522A, sc-92522B and sc-92522C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

BSPRY siRNA (h) is recommended for the inhibition of BSPRY expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

BSPRY (C-12): sc-377320 is recommended as a control antibody for monitoring of BSPRY gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor BSPRY gene expression knockdown using RT-PCR Primer: BSPRY (h)-PR: sc-92522-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.