

# Lingo-2 siRNA (h): sc-92529

## BACKGROUND

Lingo-2 (leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 2), also known as LERN3 (leucine-rich repeat neuronal protein 3) or LRRN6C, is a 606 amino acid single-pass type I membrane protein. Lingo-2 contains one immunoglobulin-like C2-type domain and 13 leucine-rich repeats. Lingo-2 is a member of the LRR gene family and has been linked to essential tremor and Parkinson's disease. Lingo-2 has a high degree of homology (61%) with Lingo-1, a regulator of neuronal death, and the expression of Lingo-2 seems to be restricted to neuronal tissue.

## REFERENCES

1. Trifunovski, A., Josephson, A., Ringman, A., Brene, S., Spenger, C. and Olson, L. 2004. Neuronal activity-induced regulation of Lingo-1. *Neuroreport* 15: 2397-2400.
2. Haines, B.P. and Rigby, P.W. 2008. Expression of the Lingo/LERN gene family during mouse embryogenesis. *Gene Expr. Patterns* 8: 79-86.
3. Benito-León, J., Louis, E.D. and Bermejo-Pareja, F. 2009. Risk of incident Parkinson's disease and parkinsonism in essential tremor: a population based study. *J. Neurol. Neurosurg. Psychiatr.* 80: 423-425.
4. Tan, E.K., Teo, Y.Y., Prakash, K.M., Li, R., Lim, H.Q., Angeles, D., Tan, L.C., Au, W.L., Yih, Y. and Zhao, Y. 2009. LINGO1 variant increases risk of familial essential tremor. *Neurology* 73: 1161-1162.
5. Clark, L.N., Park, N., Kisselev, S., Rios, E., Lee, J.H. and Louis, E.D. 2010. Replication of the LINGO1 gene association with essential tremor in a North American population. *Eur. J. Hum. Genet.* 18: 838-843.
6. Vilariño-Güell, C., Wider, C., Ross, O.A., Jasinska-Myga, B., Kachergus, J., Cobb, S.A., Soto-Ortolaza, A.I., Behrouz, B., Heckman, M.G., Diehl, N.N., Testa, C.M., Wszolek, Z.K., Uitti, R.J., et al. 2010. LINGO1 and LINGO2 variants are associated with essential tremor and Parkinson disease. *Neurogenetics* 11: 401-408.
7. Wu, Y.W., Prakash, K.M., Rong, T.Y., Li, H.H., Xiao, Q., Tan, L.C., Au, W.L., Ding, J.Q., Chen, S.D. and Tan, E.K. 2011. Lingo2 variants associated with essential tremor and Parkinson's disease. *Hum. Genet.* 129: 611-615.

## CHROMOSOMAL LOCATION

Genetic locus: LINGO2 (human) mapping to 9p21.2.

## PRODUCT

Lingo-2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Lingo-2 shRNA Plasmid (h): sc-92529-SH and Lingo-2 shRNA (h) Lentiviral Particles: sc-92529-V as alternate gene silencing products.

For independent verification of Lingo-2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-92529A, sc-92529B and sc-92529C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Lingo-2 siRNA (h) is recommended for the inhibition of Lingo-2 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Lingo-2 gene expression knockdown using RT-PCR Primer: Lingo-2 (h)-PR: sc-92529-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.