

HABP4 siRNA (h): sc-92886

BACKGROUND

The intracellular hyaluronan-binding protein 4 (HABP4 or iHABP4), also designated Ki-1/57, is a cytoplasmic and nuclear protein that is phosphorylated on serine and threonine residues. HABP4 is expressed as a protein of 413 amino acids and is phosphorylated by phorbol 12-myristate 13-acetate (PMA)-activated PKC isoforms at Thr 354 and Thr 375. HABP4 interacts with hyaluron, the adaptor protein RACK1, the transcription factor MEF2C and the chromatin remodeling factor CHD3, suggesting that it may be involved in the regulation of transcription. However, its actual function has yet to be elucidated. HABP4 is highly expressed in brain, heart, and kidney and is moderately expressed in skeletal muscle. Also, HABP4 is expressed in a variety of tumor cell lines and in activated, but not resting, leukocytes.

REFERENCES

1. Kobarg, J., Schnittger, S., Fonatsch, C., Lemke, H., Bowen, M.A., Buck, F. and Hansen, H.P. 1997. Characterization, mapping and partial cDNA sequence of the 57-kD intracellular Ki-1 antigen. *Exp. Clin. Immunogenet.* 14: 273-280.
2. Huang, L., Grammatikakis, N., Yoneda, M., Banerjee, S.D. and Toole, B.P. 2000. Molecular characterization of a novel intracellular hyaluronan-binding protein. *J. Biol. Chem.* 275: 29829-29839.
3. Lemos, T.A., Passos, D.O., Nery, F.C. and Kobarg, J. 2003. Characterization of a new family of proteins that interact with the C-terminal region of the chromatin-remodeling factor CHD-3. *FEBS Lett.* 533: 14-20.
4. Nery, F.C., Passos, D.O., Garcia, V.S. and Kobarg, J. 2004. Ki-1/57 interacts with RACK1 and is a substrate for the phosphorylation by phorbol 12-myristate 13-acetate-activated protein kinase C. *J. Biol. Chem.* 279: 11444-11455.
5. Kobarg, C.B., Kobarg, J., Crosara-Alberto, D.P., Theizen, T.H. and Franchini, K.G. 2005. MEF2C DNA-binding activity is inhibited through its interaction with the regulatory protein Ki-1/57. *FEBS Lett.* 579: 2615-2622.
6. Nery, F.C., Rui, E., Kuniyoshi, T.M. and Kobarg, J. 2006. Evidence for the interaction of the regulatory protein Ki-1/57 with p53 and its interacting proteins. *Biochem. Biophys. Res. Commun.* 341: 847-855.
7. Passos, D.O., Bressan, G.C., Nery, F.C. and Kobarg, J. 2006. Ki-1/57 interacts with PRMT1 and is a substrate for arginine methylation. *FEBS J.* 273: 3946-3961.
8. Bressan, G.C., Silva, J.C., Borges, J.C., Dos Passos, D.O., Ramos, C.H., Torriani, I.L. and Kobarg, J. 2008. Human regulatory protein Ki-1/57 has characteristics of an intrinsically unstructured protein. *J. Proteome Res.* 7: 4465-4474.
9. Bressan, G.C., Quaresma, A.J., Moraes, E.C., Manfiolli, A.O., Passos, D.O., Gomes, M.D. and Kobarg, J. 2009. Functional association of human Ki-1/57 with pre-mRNA splicing events. *FEBS J.* 276: 3770-3783.

CHROMOSOMAL LOCATION

Genetic locus: HABP4 (human) mapping to 9q22.32.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

HABP4 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HABP4 shRNA Plasmid (h): sc-92886-SH and HABP4 shRNA (h) Lentiviral Particles: sc-92886-V as alternate gene silencing products.

For independent verification of HABP4 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-92886A, sc-92886B and sc-92886C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

HABP4 siRNA (h) is recommended for the inhibition of HABP4 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor HABP4 gene expression knockdown using RT-PCR Primer: HABP4 (h)-PR: sc-92886-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.