

TBC1D3C siRNA (h): sc-93595

BACKGROUND

GTPase-activating proteins (GAPs) accelerate the intrinsic rate of GTP hydrolysis of Ras-related proteins, resulting in downregulation of their active form. TBC1D3C (TBC1 domain family, member 3C) is a 549 amino acid protein that contains one Rab-GAP TBC domain and belongs to a family of TBC domain-containing proteins that may convey a role in the regulation of cell growth and differentiation. Expressed in thymus, pancreas and testis, TBC1D3C is thought to function as a GTPase-activating protein for Rab 5, possibly participating in Rab-mediated cell cycle control and, ultimately, cellular differentiation. The gene encoding TBC1D3C maps to human chromosome 17, which comprises over 2.5% of the human genome and encodes over 1,200 genes.

REFERENCES

1. Pei, L., Peng, Y., Yang, Y., Ling, X.B., Van Eynhoven, W.G., Nguyen, K.C., Rubin, M., Hoey, T., Powers, S. and Li, J. 2002. PRC17, a novel oncogene encoding a Rab GTPase-activating protein, is amplified in prostate cancer. *Cancer Res.* 62: 5420-5424.
2. Paulding, C.A., Ruvoilo, M. and Haber, D.A. 2003. The Tre2 (USP6) oncogene is a hominoid-specific gene. *Proc. Natl. Acad. Sci. USA* 100: 2507-2511.
3. Hodzic, D., Kong, C., Wainszelbaum, M.J., Charron, A.J., Su, X. and Stahl, P.D. 2006. TBC1D3, a hominoid oncoprotein, is encoded by a cluster of paralogues located on chromosome 17q12. *Genomics* 88: 731-736.
4. Online Mendelian Inheritance in Man, OMIM[™]. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 610806. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Wainszelbaum, M.J., Charron, A.J., Kong, C., Kirkpatrick, D.S., Srikanth, P., Barbieri, M.A., Gygi, S.P. and Stahl, P.D. 2008. The hominoid-specific oncogene TBC1D3 activates Ras and modulates epidermal growth factor receptor signaling and trafficking. *J. Biol. Chem.* 283: 13233-13242.
6. Frittoli, E., Palamidessi, A., Pizzigoni, A., Lanzetti, L., Garrè, M., Troglio, F., Troilo, A., Fukuda, M., Di Fiore, P.P., Scita, G. and Confalonieri, S. 2008. The primate-specific protein TBC1D3 is required for optimal macropinocytosis in a novel ARF6-dependent pathway. *Mol. Biol. Cell* 19: 1304-1316.

CHROMOSOMAL LOCATION

Genetic locus: TBC1D3C (human) mapping to 17q12.

PRODUCT

TBC1D3C siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TBC1D3C shRNA Plasmid (h): sc-93595-SH and TBC1D3C shRNA (h) Lentiviral Particles: sc-93595-V as alternate gene silencing products.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TBC1D3C siRNA (h) is recommended for the inhibition of TBC1D3C expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

TBC1D3/3B/3C/3G (H-7): sc-376073 is recommended as a control antibody for monitoring of TBC1D3C gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TBC1D3C gene expression knockdown using RT-PCR Primer: TBC1D3C (h)-PR: sc-93595-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.