



## EVI2B siRNA (h): sc-93673

### BACKGROUND

EVI2B (ecotropic viral integration site 2B protein homolog) is a 448 amino acid protein which functions in the differentiation of melanocytes and keratinocytes. Lying within an intron of the Neurofibromin gene, the gene encoding EVI2B is transcribed from the telomere toward the centromere, which is opposite the transcription direction of the Neurofibromin gene. EVI2B is a single-pass transmembrane protein containing an extracellular domain with four glycosylation sites, a N-terminal signal peptide, a cytoplasmic hydrophilic domain and a hydrophobic transmembrane domain. Due to evidence suggesting that gene encoding the mouse homolog lies within a viral integration site that has been identified in retrovirus-induced myeloid tumors, the gene encoding EVI2B may function as an oncogene in these tumor types. With expression in peripheral blood mononuclear cells, fibroblasts, bone marrow and EBV-transformed lymphoblastoid cell lines, EVI2B is implicated in leukemogenesis.

### REFERENCES

1. Cawthon, R.M., Andersen, L.B., Buchberg, A.M., Xu, G.F., O'Connell, P., Viskochil, D., Weiss, R.B., Wallace, M.R., Marchuk, D.A. and Culver, M. 1991. cDNA sequence and genomic structure of EVI2B, a gene lying within an intron of the neurofibromatosis type 1 gene. *Genomics* 9: 446-460.
2. Viskochil, D., Cawthon, R., O'Connell, P., Xu, G.F., Stevens, J., Culver, M., Carey, J. and White, R. 1991. The gene encoding the oligodendrocyte-myelin glycoprotein is embedded within the neurofibromatosis type 1 gene. *Mol. Cell. Biol.* 11: 906-912.
3. Rasmussen, S.A., Ho, V.T., Colman, S.D., Abernathy, C.R. and Wallace, M.R. 1996. A multiplex-PCR test for EVI2A and EVI2B polymorphisms within the human NF1 gene. *Mamm. Genome* 7: 233-234.
4. Kaufmann, D., Gruener, S., Braun, F., Stark, M., Griesser, J., Hoffmeyer, S. and Bartelt, B. 1999. EVI2B, a gene lying in an intron of the neurofibromatosis type 1 (NF1) gene, is as the NF1 gene involved in differentiation of melanocytes and keratinocytes and is overexpressed in cells derived from NF1 neurofibromas. *DNA Cell Biol.* 18: 345-356.
5. Aalto, Y., El-Rifa, W., Vilpo, L., Ollila, J., Nagy, B., Vihinen, M., Vilpo, J. and Knuutila, S. 2001. Distinct gene expression profiling in chronic lymphocytic leukemia with 11q23 deletion. *Leukemia* 15: 1721-1728.
6. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 158381. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Jenne, D.E., Tinschert, S., Dorschner, M.O., Hameister, H., Stephens, K. and Kehrer-Sawatzki, H. 2003. Complete physical map and gene content of the human NF1 tumor suppressor region in human and mouse. *Genes Chromosomes Cancer* 37: 111-120.

### CHROMOSOMAL LOCATION

Genetic locus: EVI2B (human) mapping to 17q11.2.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

### PRODUCT

EVI2B siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see EVI2B shRNA Plasmid (h): sc-93673-SH and EVI2B shRNA (h) Lentiviral Particles: sc-93673-V as alternate gene silencing products.

For independent verification of EVI2B (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-93673A, sc-93673B and sc-93673C.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

EVI2B siRNA (h) is recommended for the inhibition of EVI2B expression in human cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor EVI2B gene expression knockdown using RT-PCR Primer: EVI2B (h)-PR: sc-93673-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.