

HSPB9 siRNA (h): sc-93978

BACKGROUND

The heat shock proteins (HSPs) comprise a group of highly conserved, abundantly expressed proteins with diverse functions, including the assembly and sequestering of multiprotein complexes, transportation of nascent polypeptide chains across cellular membranes and regulation of protein folding. Also referred to as molecular chaperones, heat shock proteins fall into six general families: HSP 90, HSP 70, HSP 60, the small molecular weight HSPs (like HSP 20), the immunophilins and the HSP 110 family. HSPB9 (heat shock protein β -9), also known as CT51 (cancer/testis antigen 51), is a 159 amino acid protein that belongs to the small heat shock protein (HSP20) family. Localizing to the cytoplasm as well as the nucleus, HSPB9 is expressed specifically in the testis. The gene encoding HSPB9 maps to human chromosome 17q21.2.

REFERENCES

1. Kappe, G., Verschuure, P., Philipsen, R.L., Staalduinen, A.A., Van de Boogaart, P., Boelens, W.C. and De Jong, W.W. 2001. Characterization of two novel human small heat shock proteins: protein kinase-related HSPB8 and testis-specific HSPB9. *Biochim. Biophys. Acta* 1520: 1-6.
2. de Wit, N.J., Verschuure, P., Kappe, G., King, S.M., de Jong, W.W., van Muijen, G.N. and Boelens, W.C. 2004. Testis-specific human small heat shock protein HSPB9 is a cancer/testis antigen, and potentially interacts with the dynein subunit TCTEL1. *Eur. J. Cell Biol.* 83: 337-345.
3. Heinen, R.C., Diniz-Mendes, L., Silva, J.T. and Paschoalin, V.M. 2006. Identification of the divergent calmodulin binding motif in yeast Ssb1/Hsp75 protein and in other HSP70 family members. *Braz. J. Med. Biol. Res.* 39: 1399-1408.
4. Blank, M., Bastrop, R. and Jürss, K. 2006. Stress protein response in two sibling species of Marenzelleria (Polychaeta: Spionidae): is there an influence of acclimation salinity? *Comp. Biochem. Physiol. B, Biochem. Mol. Biol.* 144: 451-462.
5. Tokalov, S.V., Pieck, S. and Gutzeit, H.O. 2007. Varying responses of human cells with discrepant p53 activity to ionizing radiation and heat shock exposure. *Cell Prolif.* 40: 24-37.
6. Im, C.N., Lee, J.S., Zheng, Y. and Seo, J.S. 2007. Iron chelation study in a normal human hepatocyte cell line suggests that tumor necrosis factor receptor-associated protein 1 (TRAP1) regulates production of reactive oxygen species. *J. Cell. Biochem.* 100: 474-486.
7. Stasyk, T., Schiefermeier, N., Skvortsov, S., Zwierzina, H., Peränen, J., Bonn, G.K. and Huber, L.A. 2007. Identification of endosomal epidermal growth factor receptor signaling targets by functional organelle proteomics. *Mol. Cell. Proteomics* 6: 908-922.
8. Vos, M.J., Kanon, B. and Kampinga, H.H. 2009. HSPB7 is a SC35 speckle resident small heat shock protein. *Biochim. Biophys. Acta* 1793: 1343-1353.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: HSPB9 (human) mapping to 17q21.2.

PRODUCT

HSPB9 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HSPB9 shRNA Plasmid (h): sc-93978-SH and HSPB9 shRNA (h) Lentiviral Particles: sc-93978-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

HSPB9 siRNA (h) is recommended for the inhibition of HSPB9 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.