



Keratin 31 siRNA (h): sc-94137

BACKGROUND

The Keratin multigene family is made of "soft" epithelial cytokeratins and "hard" hair Keratins. While the epithelial cytokeratins are involved in the layering and formation of epithelia, the hair Keratins are responsible for creating nails and hair. There are two types of Keratins: the acidic class I Keratin proteins and the basic/neutral class II Keratin proteins. As a type I hair Keratin protein, Keratin 31, also known as KRT31, HA1, HHA1 (hair Keratin, type I Ha1), K31 or KRTHA1, contains 416 amino acids and is expressed in scalp but not hairless skin. In growing hair, Keratin 31 localizes to keratinocytes of the hair cortex but is not present in medulla or inner root sheath. The gene encoding Keratin 31 maps to human chromosome 17q21.2.

REFERENCES

1. Heid, H.W., Werner, E. and Franke, W.W. 1986. The complement of native α -keratin polypeptides of hair-forming cells: a subset of eight polypeptides that differ from epithelial cytokeratins. *Differentiation* 32: 101-119.
2. Fink, P., Rogers, M.A., Korge, B., Winter, H. and Schweizer, J. 1995. A cDNA encoding the human type I hair keratin hHa1. *Biochim. Biophys. Acta* 1264: 12-14.
3. Rogers, M.A., Nischt, R., Korge, B., Krieg, T., Fink, T.M., Lichter, P., Winter, H. and Schweizer, J. 1995. Sequence data and chromosomal localization of human type I and type II hair keratin genes. *Exp. Cell Res.* 220: 357-362.
4. Rogers, M.A., Langbein, L., Praetzel, S., Moll, I., Krieg, T., Winter, H. and Schweizer, J. 1997. Sequences and differential expression of three novel human type-II hair keratins. *Differentiation* 61: 187-194.
5. Winter, H., Hofmann, I., Langbein, L., Rogers, M.A. and Schweizer, J. 1997. A splice site mutation in the gene of the human type I hair keratin hHa1 results in the expression of a tailless keratin isoform. *J. Biol. Chem.* 272: 32345-32352.
6. Bowden, P.E., Hainey, S.D., Parker, G., Jones, D.O., Zimonjic, D., Popescu, N. and Hodgins, M.B. 1998. Characterization and chromosomal localization of human hair-specific keratin genes and comparative expression during the hair growth cycle. *J. Invest. Dermatol.* 110: 158-164.
7. Langbein, L., Rogers, M.A., Winter, H., Praetzel, S., Beckhaus, U., Rackwitz, H.R. and Schweizer, J. 1999. The catalog of human hair keratins. I. Expression of the nine type I members in the hair follicle. *J. Biol. Chem.* 274: 19874-19884.
8. Online Mendelian Inheritance in Man, OMIM™. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 601077. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: KRT31 (human) mapping to 17q21.2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Keratin 31 siRNA (h) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Keratin 31 shRNA Plasmid (h): sc-94137-SH and Keratin 31 shRNA (h) Lentiviral Particles: sc-94137-V as alternate gene silencing products.

For independent verification of Keratin 31 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-94137A and sc-94137B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Keratin 31 siRNA (h) is recommended for the inhibition of Keratin 31 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Keratin 31 gene expression knockdown using RT-PCR Primer: Keratin 31 (h)-PR: sc-94137-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.