

TSPYL6 siRNA (h): sc-94760

BACKGROUND

TSPYL6 (testis-specific Y-encoded-like protein 6) is a 410 amino acid member of the nucleosome assembly protein (NAP) family. TSPYL6 is believed to be similar to Testis-specific Y-encoded protein 1 in form and function. The gene that encodes TSPYL6 is found on chromosome 2 which consists of 237 million bases encoding over 1,400 genes and making up approximately 8% of the human genome. A number of genetic diseases are linked to genes on chromosome 2. The lipid metabolic disorder sitosterolemia is associated with ABCG5 and ABCG8. An extremely rare recessive genetic disorder, Alström syndrome is due to mutations in the ALMS1 gene. Interestingly, chromosome 2 contains what appears to be a vestigial second centromere and vestigial telomeres which gives credence to the hypothesis that human chromosome 2 is the result of an ancient fusion of two ancestral chromosomes seen in modern form today in apes.

REFERENCES

1. Kido, T. and Lau, Y.F. 2006. The rat Tspy is preferentially expressed in elongated spermatids and interacts with the core histones. *Biochem. Biophys. Res. Commun.* 350: 56-67.
2. Kido, T. and Lau, Y.F. 2008. The human Y-encoded testis-specific protein interacts functionally with eukaryotic translation elongation factor eEF1A, a putative oncoprotein. *Int. J. Cancer* 123: 1573-1585.
3. Li, Y. and Lau, Y.F. 2008. TSPY and its X-encoded homologue interact with cyclin B but exert contrasting functions on cyclin-dependent kinase 1 activities. *Oncogene* 27: 6141-6150.
4. Lau, Y.F., Li, Y. and Kido, T. 2009. Gonadoblastoma locus and the TSPY gene on the human Y chromosome. *Birth Defects Res. C Embryo Today* 87: 114-122.
5. Vallangeon, B.D., Eble, J.N. and Ulbright, T.M. 2010. Macroscopic sertoli cell nodule: a study of 6 cases that presented as testicular masses. *Am. J. Surg. Pathol.* 34: 1874-1880.
6. Nickkholgh, B., Noordam, M.J., Hovingh, S.E., van Pelt, A.M., van der Veen, F. and Repping, S. 2010. Y chromosome TSPY copy numbers and semen quality. *Fertil. Steril.* 94: 1744-1747.
7. Kumar, S., Gupta, R., Kumar, S. and Ali, S. 2011. Molecular mining of alleles in water buffalo *Bubalus bubalis* and characterization of the TSPY1 and COL6A1 genes. *PLoS ONE* 6: e24958.

CHROMOSOMAL LOCATION

Genetic locus: TSPYL6 (human) mapping to 2p16.2.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

TSPYL6 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TSPYL6 shRNA Plasmid (h): sc-94760-SH and TSPYL6 shRNA (h) Lentiviral Particles: sc-94760-V as alternate gene silencing products.

For independent verification of TSPYL6 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-94760A, sc-94760B and sc-94760C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TSPYL6 siRNA (h) is recommended for the inhibition of TSPYL6 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TSPYL6 gene expression knockdown using RT-PCR Primer: TSPYL6 (h)-PR: sc-94760-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.