

# RNF103 siRNA (h): sc-94784

## BACKGROUND

RNF103 (E3 ubiquitin-protein ligase RNF103, KF-1) is a 685 amino acid multi-pass membrane protein that contains one RING-type zinc finger. RNF103 localizes to the membrane of the endoplasmic reticulum and likely interacts with DERL1 and VCP. RNF103 is believed to act as an E2-dependent E3 ubiquitin-protein ligase and is probably involved in the ER-associated protein degradation pathway. RNF103 is highly expressed in the normal cerebellum but not in the cerebral cortex. The gene encoding RNF103 maps to human chromosome 2p11.2. As the second largest human chromosome, chromosome 2 makes up approximately 8% of the human genome and contains 237 million bases encoding over 1,400 genes. A number of genetic diseases are linked to genes on chromosome 2.

## REFERENCES

1. Yasojima, K., Tsujimura, A., Mizuno, T., Shigeyoshi, Y., Inazawa, J., Kikuno, R., Kuma, K., Ohkubo, K., Hosokawa, Y., Ibata, Y., Abe, T., Miyata, T., Matsubara, K., Nakajima, K. and Hashimoto-Gotoh, T. 1997. Cloning of human and mouse cDNAs encoding novel zinc finger proteins expressed in cerebellum and hippocampus. *Biochem. Biophys. Res. Commun.* 231: 481-487.
2. Lorick, K.L., Jensen, J.P., Fang, S., Ong, A.M., Hatakeyama, S. and Weissman, A.M. 1999. RING fingers mediate ubiquitin-conjugating enzyme (E2)-dependent ubiquitination. *Proc. Natl. Acad. Sci. USA* 96: 11364-11369.
3. Yamada, M., Yamada, M., Yamazaki, S., Takahashi, K., Nishioka, G., Kudo, K., Ozawa, H., Yamada, S., Kiuchi, Y., Kamijima, K., Higuchi, T. and Momose, K. 2000. Identification of a novel gene with RING-H2 finger motif induced after chronic antidepressant treatment in rat brain. *Biochem. Biophys. Res. Commun.* 278: 150-157.
4. Nishioka, G., Yamada, M., Kudo, K., Takahashi, K., Kiuchi, Y., Higuchi, T., Momose, K., Kamijima, K. and Yamada, M. 2003. Induction of kf-1 after repeated electroconvulsive treatment and chronic antidepressant treatment in rat frontal cortex and hippocampus. *J. Neural. Transm.* 110: 277-285.
5. Hillier, L.W., Graves, T.S., Miller, W., Eichler, E.E., Bork, P., Suyama, M., Torrents, D., Waterston, R.H. and Wilson, R.K. 2005. Generation and annotation of the DNA sequences of human chromosomes 2 and 4. *Nature* 434: 724-731.
6. Maruyama, Y., Yamada, M., Takahashi, K. and Yamada, M. 2008. Ubiquitin ligase Kf-1 is involved in the endoplasmic reticulum-associated degradation pathway. *Biochem. Biophys. Res. Commun.* 374: 737-741.
7. Hashimoto-Gotoh, T., Iwabe, N., Tsujimura, A., Nakagawa, M. and Marunaka, Y. 2011. KF-1 ubiquitin ligase: anxiety suppressor model. *Cell Biochem. Biophys.* 60: 69-75.

## CHROMOSOMAL LOCATION

Genetic locus: RNF103 (human) mapping to 2p11.2.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

RNF103 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RNF103 shRNA Plasmid (h): sc-94784-SH and RNF103 shRNA (h) Lentiviral Particles: sc-94784-V as alternate gene silencing products.

For independent verification of RNF103 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-94784A, sc-94784B and sc-94784C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

RNF103 siRNA (h) is recommended for the inhibition of RNF103 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RNF103 gene expression knockdown using RT-PCR Primer: RNF103 (h)-PR: sc-94784-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.