SANTA CRUZ BIOTECHNOLOGY, INC.

Reg Iβ siRNA (h): sc-95000



BACKGROUND

The regeneration (Reg) family consists of secretory proteins involved in liver, pancreatic, gastric and intestinal cell proliferation or differentiation. Members of the REG family are divided into four subclasses, designated types I, II, III and IV. They have been implicated in the regulation of cell growth, tumorigenesis and the progression of cancer. Reg I β (regenerating islet-derived 1 β), also known as REGH, REGL or PSPS2, is a 166 amino acid pancreatic protein that contains one C-type lectin domain and belongs to subclass I of the Reg family. The gene encoding Reg I β is tandemly clustered with other Reg proteins on a region of human chromosome 2, suggesting that several Reg proteins may have arisen from a gene duplication event.

REFERENCES

- Gharib, B., Fox, M.F., Bartoli, C., Giorgi, D., Sansonetti, A., Swallow, D.M., Dagorn, J.C. and Berge-Lefranc, J.L. 1993. Human regeneration protein/ lithostathine genes map to chromosome 2p12. Ann. Hum. Genet. 57: 9-16.
- 2. Bartoli, C., Gharib, B., Giorgi, D., Sansonetti, A., Dagorn, J.C. and Berge-Lefranc, J.L. 1993. A gene homologous to the Reg gene is expressed in the human pancreas. FEBS Lett. 327: 289-293.
- Moriizumi, S., Watanabe, T., Unno, M., Nakagawara, K., Suzuki, Y., Miyashita, H., Yonekura, H. and Okamoto, H. 1994. Isolation, structural determination and expression of a novel Reg gene, human Regl β. Biochim. Biophys. Acta 1217: 199-202.
- Miyashita, H., Nakagawara, K., Mori, M., Narushima, Y., Noguchi, N., Moriizumi, S., Takasawa, S., Yonekura, H., Takeuchi, T. and Okamoto, H. 1995. Human REG family genes are tandemly ordered in a 95-kilobase region of chromosome 2p12. FEBS Lett. 377: 429-433.
- Bartoli, C., Baeza, N., Figarella, C., Pellegrini, I. and Figarella-Branger, D. 1998. Expression of peptide-23/pancreatitis-associated protein and Reg genes in human pituitary and adenomas: comparison with other fetal and adult human tissues. J. Clin. Endocrinol. Metab. 83: 4041-4046.
- 6. Sanchez, D., Figarella, C., Marchand-Pinatel, S., Bruneau, N. and Guy-Crotte, O. 2001. Preferential expression of Reg I β gene in human adult pancreas. Biochem. Biophys. Res. Commun. 284: 729-737.
- 7. Banchuin, N., Boonyasrisawat, W., Pulsawat, P., Vannasaeng, S., Deerochanawong, C., Sriussadaporn, S., Ploybutr, S., Pasurakul, T. and Yenchitsomanus, P.T. 2002. No abnormalities of Reg1 α and Reg1 β gene associated with diabetes mellitus. Diabetes Res. Clin. Pract. 55: 105-111.
- 8. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 167771. World Wide Web URL: http://www.ncbi.nlm.nih.gov/omim/

CHROMOSOMAL LOCATION

Genetic locus: REG1B (human) mapping to 2p12.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Reg I β siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Reg I β shRNA Plasmid (h): sc-95000-SH and Reg I β shRNA (h) Lentiviral Particles: sc-95000-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNAse-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNAse-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Reg I β siRNA (h) is recommended for the inhibition of Reg I β expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Reg I β gene expression knockdown using RT-PCR Primer: Reg I β (h)-PR: sc-95000-PR (20 μ I). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.