

ALDH8A1 siRNA (h): sc-95150

BACKGROUND

Aldehyde dehydrogenases (ALDHs) mediate the NADP⁺-dependent oxidation of aldehydes into acids and play an important role in the detoxification of alcohol-derived acetaldehyde, as well as in lipid peroxidation and in the metabolism of corticosteroids, biogenic amines and neurotransmitters. ALDH8A1 (aldehyde dehydrogenase 8 family, member A1), also known as ALDH12, is a 487 amino acid protein that localizes to the cytoplasm and belongs to the aldehyde dehydrogenase family. Expressed in kidney and liver, ALDH8A1 converts 9-*cis*-retinal to 9-*cis*-retinoic acid. 9-*cis*-retinoic acid activates retinoid X receptors, a family of nuclear receptors which are involved in regulating multiple signaling pathways. Three isoforms exist due to alternative splicing events.

REFERENCES

- Heyman, R.A., Mangelsdorf, D.J., Dyck, J.A., Stein, R.B., Eichele, G., Evans, R.M. and Thaller, C. 1992. 9-*cis*-retinoic acid is a high affinity ligand for the retinoid X receptor. *Cell* 68: 397-406.
- Lin, M. and Napoli, J.L. 2000. cDNA cloning and expression of a human aldehyde dehydrogenase (ALDH) active with 9-*cis*-retinal and identification of a rat ortholog, ALDH12. *J. Biol. Chem.* 275: 40106-40112.
- Vasilou, V. and Pappa, A. 2000. Polymorphisms of human aldehyde dehydrogenases. Consequences for drug metabolism and disease. *Pharmacology* 61: 192-198.
- Zhuang, R., Lin, M. and Napoli, J.L. 2002. *cis*-Retinol/androgen dehydrogenase, isozyme 3 (CRAD3): a short-chain dehydrogenase active in a reconstituted path of 9-*cis*-retinoic acid biosynthesis in intact cells. *Biochemistry* 41: 3477-3483.
- Ahuja, H.S., Szanto, A., Nagy, L. and Davies, P.J. 2003. The retinoid X receptor and its ligands: versatile regulators of metabolic function, cell differentiation and cell death. *J. Biol. Regul. Homeost. Agents* 17: 29-45.

CHROMOSOMAL LOCATION

Genetic locus: ALDH8A1 (human) mapping to 6q23.3.

PRODUCT

ALDH8A1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ALDH8A1 shRNA Plasmid (h): sc-95150-SH and ALDH8A1 shRNA (h) Lentiviral Particles: sc-95150-V as alternate gene silencing products.

For independent verification of ALDH8A1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-95150A, sc-95150B and sc-95150C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ALDH8A1 siRNA (h) is recommended for the inhibition of ALDH8A1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

ALDH8A1 (E-2): sc-515006 is recommended as a control antibody for monitoring of ALDH8A1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ALDH8A1 gene expression knockdown using RT-PCR Primer: ALDH8A1 (h)-PR: sc-95150-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.