

connexin 25 siRNA (h): sc-95275

BACKGROUND

The connexin family of proteins form hexameric complexes called "connexons" that facilitate movement of low molecular weight proteins between cells via gap junctions. Connexin proteins share a common topology of four transmembrane α -helical domains, two extracellular loops, a cytoplasmic loop and cytoplasmic N- and C-termini. Many of the key functional differences arise from specific amino acid substitutions in the most highly conserved domains, the transmembrane and extracellular regions. Each of the approximately 20 connexin isoforms produces channels with distinct permeabilities and electrical and chemical sensitivities; therefore, one connexin usually cannot fully substitute for another. Consequently, a wide variety of malignant phenotypes associate with decreased connexin expression and gap junction communication, dependent on the particular connexin that is affected. Connexin 25, also known as GJB7 (gap junction protein, β 7) or CX25, is a 223 amino acid multi-pass membrane protein that is weakly expressed in placenta. Belonging to the connexin family and the β -type (group I) subfamily, connexin 25 is encoded by a gene located on human chromosome 6q14.3.

REFERENCES

1. von Maltzahn, J., Euwens, C., Willecke, K. and Söhl, G. 2004. The novel mouse connexin 39 gene is expressed in developing striated muscle fibers. *J. Cell Sci.* 117: 5381-5392.
2. Xia, C.H., Liu, H., Cheung, D., Cheng, C., Wang, E., Du, X., Beutler, B., Lo, W.K. and Gong, X. 2006. Diverse gap junctions modulate distinct mechanisms for fiber cell formation during lens development and cataractogenesis. *Development* 133: 2033-2040.
3. Dunia, I., Cibert, C., Gong, X., Xia, C.H., Recouvreur, M., Levy, E., Kumar, N., Bloemendal, H. and Benedetti, E.L. 2006. Structural and immunocytochemical alterations in eye lens fiber cells from Cx46 and Cx50 knockout mice. *Eur. J. Cell Biol.* 85: 729-752.
4. Xia, C.H., Cheng, C., Huang, Q., Cheung, D., Li, L., Dunia, I., Benedetti, L.E., Horwitz, J. and Gong, X. 2006. Absence of α 3 (Cx46) and α 8 (Cx50) connexins leads to cataracts by affecting lens inner fiber cells. *Exp. Eye Res.* 83: 688-696.
5. Xia, C.H., Cheung, D., DeRosa, A.M., Chang, B., Lo, W.K., White, T.W. and Gong, X. 2006. Knock-in of α 3 connexin prevents severe cataracts caused by an α 8 point mutation. *J. Cell Sci.* 119: 2138-2144.
6. Tang, Y., Liu, X., Zoltoski, R.K., Novak, L.A., Herrera, R.A., Richard, I., Kuszak, J.R. and Kumar, N.M. 2007. Age-related cataracts in α 3Cx46-knockout mice are dependent on a calpain 3 isoform. *Invest. Ophthalmol. Vis. Sci.* 48: 2685-2694.
7. Cheng, C., Xia, C.H., Li, L., White, T.W., Niimi, J. and Gong, X. 2008. Gap junction communication influences intercellular protein distribution in the lens. *Exp. Eye Res.* 86: 966-974.
8. Hoang, Q.V., Qian, H. and Ripps, H. 2010. Functional analysis of hemichannels and gap-junctional channels formed by connexins 43 and 46. *Mol. Vis.* 16: 1343-1352.

CHROMOSOMAL LOCATION

Genetic locus: GJB7 (human) mapping to 6q14.3.

PRODUCT

connexin 25 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see connexin 25 shRNA Plasmid (h): sc-95275-SH and connexin 25 shRNA (h) Lentiviral Particles: sc-95275-V as alternate gene silencing products.

For independent verification of connexin 25 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-95275A, sc-95275B and sc-95275C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

connexin 25 siRNA (h) is recommended for the inhibition of connexin 25 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor connexin 25 gene expression knockdown using RT-PCR Primer: connexin 25 (h)-PR: sc-95275-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.