

VPS37B siRNA (h): sc-95762

BACKGROUND

Vacuolar sorting proteins (VPSs) are required for proper trafficking of endocytic and biosynthetic proteins to the vacuole and play an important role in the budding process of cells. VPS37B (vacuolar protein sorting 37 homolog B), also known as ESCRT-I complex subunit VPS37B, is a 285 amino acid peripheral membrane protein and is a member of the VPS37 family. Expressed in macrophages and lymphocytes, VPS37B is recruited to the endosomal membrane in a VPS4A-dependent fashion. VPS37B is a component of the ESCRT-I complex, which is known to regulate vesicular trafficking. VPS37B is required for the sorting of endocytic ubiquitinated cargo into multivesicular bodies and may be involved in cell growth and differentiation.

REFERENCES

1. Stuchell, M.D., Garrus, J.E., Müller, B., Stray, K.M., Ghaffarian, S., McKinnon, R., Kräusslich, H.G., Morham, S.G. and Sundquist, W.I. 2004. The human endosomal sorting complex required for transport (ESCRT-I) and its role in HIV-1 budding. *J. Biol. Chem.* 279: 36059-36071.
2. Bache, K.G., Slagsvold, T., Cabezas, A., Rosendal, K.R., Raiborg, C. and Stenmark, H. 2004. The growth-regulatory protein HCRP1/hVps37A is a subunit of mammalian ESCRT-I and mediates receptor down-regulation. *Mol. Biol. Cell* 15: 4337-4346.
3. Eastman, S.W., Martin-Serrano, J., Chung, W., Zang, T. and Bieniasz, P.D. 2005. Identification of human VPS37C, a component of endosomal sorting complex required for transport-I important for viral budding. *J. Biol. Chem.* 280: 628-636.
4. Gill, D.J., Teo, H., Sun, J., Perisic, O., Veprintsev, D.B., Vallis, Y., Emr, S.D. and Williams, R.L. 2007. Structural studies of phosphoinositide 3-kinase-dependent traffic to multivesicular bodies. *Biochem. Soc. Symp.* 74: 47-57.
5. Morita, E., Sandrin, V., Alam, S.L., Eckert, D.M., Gygi, S.P. and Sundquist, W.I. 2007. Identification of human MVB12 proteins as ESCRT-I subunits that function in HIV budding. *Cell Host Microbe* 2: 41-53.
6. Morita, E., Sandrin, V., Chung, H.Y., Morham, S.G., Gygi, S.P., Rodesch, C.K. and Sundquist, W.I. 2007. Human ESCRT and ALIX proteins interact with proteins of the midbody and function in cytokinesis. *EMBO J.* 26: 4215-4227.

CHROMOSOMAL LOCATION

Genetic locus: VPS37B (human) mapping to 12q24.31.

PRODUCT

VPS37B siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see VPS37B shRNA Plasmid (h): sc-95762-SH and VPS37B shRNA (h) Lentiviral Particles: sc-95762-V as alternate gene silencing products.

For independent verification of VPS37B (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-95762A, sc-95762B and sc-95762C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

VPS37B siRNA (h) is recommended for the inhibition of VPS37B expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

VPS37B (D-4): sc-390144 is recommended as a control antibody for monitoring of VPS37B gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor VPS37B gene expression knockdown using RT-PCR Primer: VPS37B (h)-PR: sc-95762-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.