

OY-TES-1 siRNA (h): sc-95804

BACKGROUND

OY-TES-1, also known as ACRBP (Acrosin-binding protein) or SP32, is a 543 amino acid secreted protein that is found on sperm acrosomes, where it colocalizes with Acrosin. Expressed in normal testicular tissue, OY-TES-1 binds to Acrosin and is thought to mediate the packaging and condensation of Acrosin in the acrosomal matrix. In addition to its expression in testis, OY-TES-1 is present in cancer tissue throughout the body, including breast, bladder, liver and lung carcinomas, suggesting an involvement in tumor formation and metastasis. The gene encoding OY-TES-1 maps to human chromosome 12, which encodes over 1,100 genes and comprises approximately 4.5% of the human genome. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

1. Baba, T., et al. 1994. An acrosomal protein, sp32, in mammalian sperm is a binding protein specific for two proacrosins and an Acrosin intermediate. *J. Biol. Chem.* 269: 10133-10140.
2. Ono, T., et al. 2001. Identification of proacrosin binding protein sp32 precursor as a human cancer/testis antigen. *Proc. Natl. Acad. Sci. USA* 98: 3282-3287.
3. Tümer, Z., et al. 2002. Genomic structure, chromosome mapping and expression analysis of the human AVIL gene, and its exclusion as a candidate for locus for inflammatory bowel disease at 12q13-14 (IBD2). *Gene* 288: 179-185.
4. Dube, C., et al. 2005. The proacrosin binding protein, sp32, is tyrosine phosphorylated during capacitation of pig sperm. *J. Androl.* 26: 519-528.
5. Okumura, H., et al. 2005. Identification of an HLA-A24-restricted OY-TES-1 epitope recognized by cytotoxic T-cells. *Microbiol. Immunol.* 49: 1009-1016.
6. Online Mendelian Inheritance in Man, OMIM[™]. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 608352. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. Tammela, J., et al. 2006. OY-TES-1 expression and serum immunoreactivity in epithelial ovarian cancer. *Int. J. Oncol.* 29: 903-910.

CHROMOSOMAL LOCATION

Genetic locus: ACRBP (human) mapping to 12p13.31.

PRODUCT

OY-TES-1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see OY-TES-1 shRNA Plasmid (h): sc-95804-SH and OY-TES-1 shRNA (h) Lentiviral Particles: sc-95804-V as alternate gene silencing products.

For independent verification of OY-TES-1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-95804A, sc-95804B and sc-95804C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

OY-TES-1 siRNA (h) is recommended for the inhibition of OY-TES-1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

OY-TES-1 (G-5): sc-390594 is recommended as a control antibody for monitoring of OY-TES-1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor OY-TES-1 gene expression knockdown using RT-PCR Primer: OY-TES-1 (h)-PR: sc-95804-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.