

AMN1 siRNA (h): sc-95825

BACKGROUND

AMN1 (antagonist of mitotic exit network 1 homolog) is a 258 amino acid protein belonging to the AMN1 family. Encoded by a gene that maps to human chromosome 12p11.21, AMN1 contains twelve degenerate leucine-rich repeat (LRR) motifs and modulates Ste12 binding to promoters of multiple genes. As a daughter cell-specific protein, AMN1 inhibits ELL function to deactivate the mitotic exit state and reset the cell cycle in G₁. AMN1 does not halt mitotic exit, but turns it off afterward. Induced by ELL after exit activation, AMN1 disrupts TEM1-Cdc15 interaction by competing with Cdc15 for TEM1 binding and by facilitating inactivation of Cdc14. Defects in AMN1 can result in disruption of both the spindle assembly and nuclear orientation checkpoints. AMN1 may also mediate gene activity and phenotypic diversity.

REFERENCES

1. Wang, Y., Shirogane, T., Liu, D., Harper, J.W. and Elledge, S.J. 2003. Exit from exit: resetting the cell cycle through Amn1 inhibition of G protein signaling. *Cell* 112: 697-709.
2. Yvert, G., Brem, R.B., Whittle, J., Akey, J.M., Foss, E., Smith, E.N., Mackelprang, R. and Kruglyak, L. 2003. Trans-acting regulatory variation in *Saccharomyces cerevisiae* and the role of transcription factors. *Nat. Genet.* 35: 57-64.
3. Jansen, R.C. and Nap, J.P. 2004. Regulating gene expression: surprises still in store. *Trends Genet.* 20: 223-225.
4. Gresham, D., Ruderfer, D.M., Pratt, S.C., Schacherer, J., Dunham, M.J., Botstein, D. and Kruglyak, L. 2006. Genome-wide detection of polymorphisms at nucleotide resolution with a single DNA microarray. *Science* 311: 1932-1936.
5. Sun, W., Yu, T. and Li, K.C. 2007. Detection of eQTL modules mediated by activity levels of transcription factors. *Bioinformatics* 23: 2290-2297.
6. Zheng, W., Zhao, H., Mancera, E., Steinmetz, L.M. and Snyder, M. 2010. Genetic analysis of variation in transcription factor binding in yeast. *Nature* 464: 1187-1191.
7. Zhang, W., Zhu, J., Schadt, E.E. and Liu, J.S. 2010. A Bayesian partition method for detecting pleiotropic and epistatic eQTL modules. *PLoS Comput. Biol.* 6: e1000642.

CHROMOSOMAL LOCATION

Genetic locus: AMN1 (human) mapping to 12p11.21.

PRODUCT

AMN1 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see AMN1 shRNA Plasmid (h): sc-95825-SH and AMN1 shRNA (h) Lentiviral Particles: sc-95825-V as alternate gene silencing products.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

AMN1 siRNA (h) is recommended for the inhibition of AMN1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor AMN1 gene expression knockdown using RT-PCR Primer: AMN1 (h)-PR: sc-95825-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.