

EF-CAB4B siRNA (h): sc-96150

BACKGROUND

EF-CAB4B (EF-hand calcium-binding domain-containing protein 4B), also known as CRACR2A (calcium release-activated calcium channel regulator 2A), is a 395 amino acid protein belonging to the EF-CAB4 family. Localizing to cytoplasm, EF-CAB4B contains two EF-hand domains and exists as two alternatively spliced isoforms. At low Ca^{2+} concentrations, EF-CAB4B acts as a calcium-sensor, facilitating the clustering of Orai1 and Stim1 at the junctional regions between plasma membrane and endoplasmic reticulum, leading to regulation of CRAC channel activation. The gene encoding EF-CAB4B maps to human chromosome 12p13.32. Encoding over 1,100 genes within 132 million bases, chromosome 12 makes up about 4.5% of the human genome. A number of skeletal deformities are linked to chromosome 12 including hypochondrogenesis, achondrogenesis and Kniest dysplasia.

REFERENCES

- Allen, T.L., et al. 1996. Cytogenetic and molecular analysis in trisomy 12p. *Am. J. Med. Genet.* 63: 250-256.
- Körkkö, J., et al. 2000. Widely distributed mutations in the COL2A1 gene produce achondrogenesis type II/hypochondrogenesis. *Am. J. Med. Genet.* 92: 95-100.
- Kelley, J., et al. 2005. Comparative genomics of natural killer cell receptor gene clusters. *PLoS Genet.* 1: 129-139.
- Lim, J., et al. 2006. A protein-protein interaction network for human inherited ataxias and disorders of Purkinje cell degeneration. *Cell* 125: 801-814.
- Dwyer, E., et al. 2010. Genotype-phenotype correlation in DTDST dysplasias: Atelosteogenesis type II and diastrophic dysplasia variant in one family. *Am. J. Med. Genet. A* 152A: 3043-3050.
- Chalasani, N., et al. 2010. Genome-wide association study identifies variants associated with histologic features of nonalcoholic fatty liver disease. *Gastroenterology* 139: 1567-1576.
- Srikanth, S., et al. 2010. A novel EF-hand protein, CRACR2A, is a cytosolic Ca^{2+} sensor that stabilizes CRAC channels in T cells. *Nat. Cell Biol.* 12: 436-446.
- Jongmans, M.C., et al. 2011. Cancer risk in patients with Noonan syndrome carrying a PTPN11 mutation. *Eur. J. Hum. Genet.* 19: 870-874.

CHROMOSOMAL LOCATION

Genetic locus: EFCAB4B (human) mapping to 12p13.32.

PRODUCT

EF-CAB4B siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see EF-CAB4B shRNA Plasmid (h): sc-96150-SH and EF-CAB4B shRNA (h) Lentiviral Particles: sc-96150-V as alternate gene silencing products.

For independent verification of EF-CAB4B (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-96150A, sc-96150B and sc-96150C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

EF-CAB4B siRNA (h) is recommended for the inhibition of EF-CAB4B expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor EF-CAB4B gene expression knockdown using RT-PCR Primer: EF-CAB4B (h)-PR: sc-96150-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.