

β-1,4-GalNAc-T3 siRNA (h): sc-96247

BACKGROUND

β-1,4-GalNAc-T3 (β-1,4-N-acetyl-galactosaminyl transferase 3), also known as B4GALNT3, is a 998 amino acid single-pass type II membrane protein belonging to the chondroitin N-acetylgalactosaminyltransferase family. β-1,4-GalNAc-T3 is highly expressed in testis, colon and stomach, and weakly expressed in other tissues. β-1,4-GalNAc-T3 exhibits subcellular localization to apical Golgi and exists as two alternatively spliced isoforms. β-1,4-GalNAc-T3 is involved in the mediation of N,N'-diacetyllactosamine formation on gastric mucosa and in N-acetyl-β-glucosaminyl-glycoprotein 4-β-N-acetylgalactosaminyltransferase activities. β-1,4-GalNAc-T3 is integral to cellular membranes and may function in a protective capacity against sudden cardiac arrest.

REFERENCES

1. Sato, T., et al. 2003. Molecular cloning and characterization of a novel human β1,4-N-acetylgalactosaminyltransferase, β4GalNAc-T3, responsible for the synthesis of N,N'-diacetyllactosamine, galNAc β1-4GlcNAc. *J. Biol. Chem.* 278: 47534-47544.
2. Ikehara, Y., et al. 2006. Apical Golgi localization of N,N'-diacetyllactosamine synthase, β4GalNAc-T3, is responsible for LacdiNAc expression on gastric mucosa. *Glycobiology* 16: 777-785.
3. Nguyen, S.T., et al. 2007. Identification of a predictive gene expression signature of cervical lymph node metastasis in oral squamous cell carcinoma. *Cancer Sci.* 98: 740-746.
4. Huang, J., et al. 2007. β1,4-N-acetylgalactosaminyltransferase III enhances malignant phenotypes of colon cancer cells. *Mol. Cancer Res.* 5: 543-552.
5. Rooryck, C., et al. 2009. 2.3 Mb terminal deletion in 12p13.33 associated with oculaurovertebral spectrum and evaluation of WNT5B as a candidate gene. *Eur. J. Med. Genet.* 52: 446-449.
6. Ito, H., et al. 2009. Strategy for glycoproteomics: identification of glyco-alteration using multiple glycan profiling tools. *J. Proteome Res.* 8: 1358-1367.
7. Narimatsu, H., et al. 2010. A strategy for discovery of cancer glyco-biomarkers in serum using newly developed technologies for glycoproteomics. *FEBS J.* 277: 95-105.
8. Fukushima, K., et al. 2010. α1,2-fucosylated and β-N-acetylgalactosaminylated prostate-specific antigen as an efficient marker of prostatic cancer. *Glycobiology* 20: 452-460.
9. Arking, D.E., et al. 2010. Genome-wide association study identifies GPC5 as a novel genetic locus protective against sudden cardiac arrest. *PLoS ONE* 5: e9879.

CHROMOSOMAL LOCATION

Genetic locus: B4GALNT3 (human) mapping to 12p13.33.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

β-1,4-GalNAc-T3 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see β-1,4-GalNAc-T3 shRNA Plasmid (h): sc-96247-SH and β-1,4-GalNAc-T3 shRNA (h) Lentiviral Particles: sc-96247-V as alternate gene silencing products.

For independent verification of β-1,4-GalNAc-T3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-96247A, sc-96247B and sc-96247C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

β-1,4-GalNAc-T3 siRNA (h) is recommended for the inhibition of β-1,4-GalNAc-T3 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor β-1,4-GalNAc-T3 gene expression knockdown using RT-PCR Primer: β-1,4-GalNAc-T3 (h)-PR: sc-96247-PR (20 μl). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.