



ACRV1 siRNA (h): sc-96294

BACKGROUND

ACRV1 (acrosomal vesicle protein 1), also known as acrosomal protein SP-10 or SPACA2, is a 265 amino acid protein. ACRV1 is encoded by a gene that maps to human chromosome 11q24.2, at the junction between 11q23 and 11q24. Containing four exons, ACRV1 may experience cryptic splicing and exon skipping. ACRV1 exists as 11 alternatively spliced isoforms and may be involved in sperm-zona binding or penetration. ACRV1 encodes a testis-specific, differentiation antigen, acrosomal vesicle protein 1 that originates in the acrosomal vesicle during spermatogenesis, and is affiliated with acrosomal membranes and mature sperm matrix. ACRV1 is a potential contraceptive vaccine immunogen.

REFERENCES

1. Wright, R.M., John, E., Klotz, K., Flickinger, C.J. and Herr, J.C. 1990. Cloning and sequencing of cDNAs coding for the human intra-acrosomal antigen SP-10. *Biol. Reprod.* 42: 693-701.
2. Homyk, M., Anderson, D.J., Wolff, H. and Herr, J.C. 1990. Differential diagnosis of immature germ cells in semen utilizing monoclonal antibody MHS-10 to the intra-acrosomal antigen SP-10. *Fertil. Steril.* 53: 323-330.
3. Herr, J.C., Wright, R.M., Flickinger, C.J., Eddy, R.L. and Shows, T.B. 1991. Assignment of the gene for human intra-acrosomal protein SP-10 to the p12→q13 region of chromosome 11. *J. Androl.* 12: 281-287.
4. Herr, J.C., Klotz, K., Shannon, J., Wright, R.M. and Flickinger, C.J. 1992. Purification and microsequencing of the intra-acrosomal protein SP-10. Evidence that SP-10 heterogeneity results from endoproteolytic processes. *Biol. Reprod.* 47: 11-20.
5. Wright, R.M., Suri, A.K., Kornreich, B., Flickinger, C.J. and Herr, J.C. 1993. Cloning and characterization of the gene coding for the human acrosomal protein SP-10. *Biol. Reprod.* 49: 316-325.
6. Golden, W.L., von Kap-Herr, C., Kurth, B., Wright, R.M., Flickinger, C.J., Eddy, R., Shows, T. and Herr, J.C. 1993. Refinement of the localization of the gene for human intraacrosomal protein SP-10 (ACRV1) to the junction of bands q23→q24 of chromosome 11 by nonisotopic *in situ* hybridization. *Genomics* 18: 446-449.
7. Foster, J.A., Klotz, K.L., Flickinger, C.J., Thomas, T.S., Wright, R.M., Castillo, J.R. and Herr, J.C. 1994. Human SP-10: acrosomal distribution, processing, and fate after the acrosome reaction. *Biol. Reprod.* 51: 1222-1231.
8. Freemerman, A.J., Flickinger, C.J. and Herr, J.C. 1995. Characterization of alternatively spliced human SP-10 mRNAs. *Mol. Reprod. Dev.* 41: 100-108.
9. Reddi, P.P., Shore, A.N., Acharya, K.K. and Herr, J.C. 2002. Transcriptional regulation of spermiogenesis: insights from the study of the gene encoding the acrosomal protein SP-10. *J. Reprod. Immunol.* 53: 25-36.

CHROMOSOMAL LOCATION

Genetic locus: ACRV1 (human) mapping to 11q24.2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

ACRV1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ACRV1 shRNA Plasmid (h): sc-96294-SH and ACRV1 shRNA (h) Lentiviral Particles: sc-96294-V as alternate gene silencing products.

For independent verification of ACRV1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-96294A, sc-96294B and sc-96294C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ACRV1 siRNA (h) is recommended for the inhibition of ACRV1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ACRV1 gene expression knockdown using RT-PCR Primer: ACRV1 (h)-PR: sc-96294-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.