

Synaptotagmin IX siRNA (h): sc-96623

BACKGROUND

Synaptotagmins are a large family of synaptic vesicle type III integral membrane proteins that function as regulators of both exocytosis and endocytosis and are involved in neurotransmitter secretion from small secretory vesicles. Synaptotagmin IX, also known as SYT9 (Synaptotagmin-9), is a 491 amino acid protein that localizes to the membrane. Like other Synaptotagmin proteins, Synaptotagmin IX is involved in the calcium-dependent exocytosis of secretory vesicles and is thought to act as a calcium sensor during vesicular trafficking. Synaptotagmin IX contains two C2 domains through which it can bind three calcium ions per subunit. It has been suggested that Synaptotagmin IX is required for the Ca^{++} -dependent release of norepinephrine.

REFERENCES

1. Perin, M.S. 1996. Mirror image motifs mediate the interaction of the COOH terminus of multiple synaptotagmins with the neurexins and calmodulin. *Biochemistry* 35: 13808-13816.
2. Mizutani, A., Fukuda, M., Ibata, K., Shiraishi, Y. and Mikoshiba, K. 2000. SYNCRIP, a cytoplasmic counterpart of heterogeneous nuclear ribonucleoprotein R, interacts with ubiquitous synaptotagmin isoforms. *J. Biol. Chem.* 275: 9823-9831.
3. Dubois, T., Howell, S., Zemlickova, E. and Aitken, A. 2002. Identification of casein kinase α interacting protein partners. *FEBS Lett.* 517: 167-171.
4. Fukuda, M., Kowalchuk, J.A., Zhang, X., Martin, T.F. and Mikoshiba, K. 2002. Synaptotagmin IX regulates Ca^{++} -dependent secretion in PC-12 cells. *J. Biol. Chem.* 277: 4601-4604.
5. Zhang, X., Kim-Miller, M.J., Fukuda, M., Kowalchuk, J.A. and Martin, T.F. 2002. Ca^{++} -dependent synaptotagmin binding to SNAP 25 is essential for Ca^{++} -triggered exocytosis. *Neuron* 34: 599-611.
6. Haberman, Y., Grimberg, E., Fukuda, M. and Sagi-Eisenberg, R. 2003. Synaptotagmin IX, a possible linker between the perinuclear endocytic recycling compartment and the microtubules. *J. Cell Sci.* 116: 4307-4318.
7. Morenilla-Palao, C., Planells-Cases, R., García-Sanz, N. and Ferrer-Montiel, A. 2004. Regulated exocytosis contributes to protein kinase C potentiation of vanilloid receptor activity. *J. Biol. Chem.* 279: 25665-25672.

CHROMOSOMAL LOCATION

Genetic locus: SYT9 (human) mapping to 11p15.4.

PRODUCT

Synaptotagmin IX siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Synaptotagmin IX shRNA Plasmid (h): sc-96623-SH and Synaptotagmin IX shRNA (h) Lentiviral Particles: sc-96623-V as alternate gene silencing products.

For independent verification of Synaptotagmin IX (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-96623A, sc-96623B and sc-96623C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Synaptotagmin IX siRNA (h) is recommended for the inhibition of Synaptotagmin IX expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Synaptotagmin IX (E-11): sc-398592 is recommended as a control antibody for monitoring of Synaptotagmin IX gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Synaptotagmin IX gene expression knockdown using RT-PCR Primer: Synaptotagmin IX (h)-PR: sc-96623-PR (20 μ l). Annealing temperature for the primers should be $55-60^{\circ}$ C and the extension temperature should be $68-72^{\circ}$ C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.