

# HRASLS2 siRNA (h): sc-96654

## BACKGROUND

Phospholipase A<sub>2</sub>s (PLA<sub>2</sub>s) constitute a family of esterases that hydrolyze the sn-2-acyl ester bond in glycerophospholipid molecules. These enzymes are generally calcium-dependent and have been found both intra- and extracellularly. By hydrolyzing the sn-2 bond in glycerophospholipids, PLA<sub>2</sub>s release fatty acids. One such fatty acid, arachidonic acid, generates substrates for the initiation of the arachidonic acid cascade that produces various eicosanoids, many of which are potent mediators of inflammation. HRASLS2 (HRAS-like suppressor 2), also known as PLA1/2-2, is a 162 amino acid cytoplasmic protein that belongs to the H-rev107 family. Expressed in liver, kidney, small intestine and colon, HRASLS2 exhibits PLA1/2 activity by catalyzing the calcium-independent hydrolysis of acyl groups in various phosphatidylcholines (PCs) and phosphatidylethanolamine (PE).

## REFERENCES

1. Johnson, L.K., Frank, S., Vades, P., Pruzanski, W., Lusis, A.J. and Seilhamer, J.J. 1990. Localization and evolution of two human phospholipase A<sub>2</sub> genes and two related genetic elements. *Adv. Exp. Med. Biol.* 275: 17-34.
2. Tischfield, J.A., Xia, Y.R., Shih, D.M., Klisak, I., Chen, J., Engle, S.J., Siakotos, A.N., Winstead, M.V., Seilhamer, J.J., Allamand, V., Gyapay, G. and Lusis, A.J. 1996. Low-molecular-weight, calcium-dependent phospholipase A<sub>2</sub> genes are linked and map to homologous chromosome regions in mouse and human. *Genomics* 32: 328-333.
3. Lambeau, G. and Lazdunski, M. 1999. Receptors for a growing family of secreted phospholipases A<sub>2</sub>. *Trends Pharmacol. Sci.* 20: 162-170.
4. Shyu, R.Y., Hsieh, Y.C., Tsai, F.M., Wu, C.C. and Jiang, S.Y. 2008. Cloning and functional characterization of the HRASLS2 gene. *Amino Acids* 35: 129-137.
5. Uyama, T., Jin, X.H., Tsuboi, K., Tonai, T. and Ueda, N. 2009. Characterization of the human tumor suppressors TIG3 and HRASLS2 as phospholipid-metabolizing enzymes. *Biochim. Biophys. Acta* 1791: 1114-1124.
6. Abbott, M.J., Tang, T. and Sul, H.S. 2010. The Role of Phospholipase A<sub>2</sub>-derived Mediators in Obesity. *Drug Discov Today Dis Mech.* 7: e213-e218.

## CHROMOSOMAL LOCATION

Genetic locus: PLAAT2 (human) mapping to 11q12.3.

## PRODUCT

HRASLS2 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HRASLS2 shRNA Plasmid (h): sc-96654-SH and HRASLS2 shRNA (h) Lentiviral Particles: sc-96654-V as alternate gene silencing products.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

HRASLS2 siRNA (h) is recommended for the inhibition of HRASLS2 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor HRASLS2 gene expression knockdown using RT-PCR Primer: HRASLS2 (h)-PR: sc-96654-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.