

## ▶ TRIM48 siRNA (h): sc-96902

### BACKGROUND

The tripartite motif (TRIM) family of proteins are characterized by a conserved TRIM domain that includes a coiled-coil region, a B-box type zinc finger, one RING finger and three zinc-binding domains. Many TRIM proteins are induced by type I and type II interferons, making them crucial for development of pathogen-resistance. TRIM48 (tripartite motif-containing protein 48), also known as RING finger protein 101, is a 208 amino acid protein that contains one RING-type zinc finger and one B box-type zinc finger. The gene encoding TRIM48 maps to human chromosome 11, which houses over 1,400 genes and comprises nearly 4% of the human genome. Jervell and Lange-Nielsen syndrome, Jacobsen syndrome, Niemann-Pick disease, hereditary angioedema and Smith-Lemli-Opitz syndrome are associated with defects in genes that maps to chromosome 11.

### REFERENCES

1. Meroni, G. and Diez-Roux, G. 2005. TRIM/RBCC, a novel class of "single protein RING finger" E3 ubiquitin ligases. *Bioessays* 27: 1147-1157.
2. Nisole, S., Stoye, J.P. and Saïb, A. 2005. TRIM family proteins: retroviral restriction and antiviral defence. *Nat. Rev. Microbiol.* 3: 799-808.
3. Towers, G.J. 2007. The control of viral infection by tripartite motif proteins and cyclophilin A. *Retrovirology* 4: 40.
4. Ozato, K., Shin, D.M., Chang, T.H. and Morse, H.C. 2008. TRIM family proteins and their emerging roles in innate immunity. *Nat. Rev. Immunol.* 8: 849-860.
5. Du Pasquier, L. 2009. Fish "n" TRIMs. *J. Biol.* 8: 50.
6. Yang, K., Shi, H.X., Liu, X.Y., Shan, Y.F., Wei, B., Chen, S. and Wang, C. 2009. TRIM21 is essential to sustain IFN regulatory factor 3 activation during antiviral response. *J. Immunol.* 182: 3782-3792.
7. Carthagena, L., Bergamaschi, A., Luna, J.M., David, A., Uchil, P.D., Margottin-Goguet, F., Mothes, W., Hazan, U., Transy, C., Pancino, G. and Nisole, S. 2009. Human TRIM gene expression in response to interferons. *PLoS ONE* 4: e4894.
8. Munir, M. 2010. TRIM proteins: another class of viral victims. *Sci. Signal.* 3: jc2.

### CHROMOSOMAL LOCATION

Genetic locus: TRIM48 (human) mapping to 11q11.

### PRODUCT

TRIM48 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TRIM48 shRNA Plasmid (h): sc-96902-SH and TRIM48 shRNA (h) Lentiviral Particles: sc-96902-V as alternate gene silencing products.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

TRIM48 siRNA (h) is recommended for the inhibition of TRIM48 expression in human cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TRIM48 gene expression knockdown using RT-PCR Primer: TRIM48 (h)-PR: sc-96902-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.