



ZNF45 siRNA (h): sc-97422

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which is thought to interact with KAP1, thereby recruiting histone modifying proteins. As a member of the krueppel C₂H₂-type zinc-finger protein family, ZNF45 (zinc finger protein 45), also known as BRC1744, Zinc finger protein 13 and zinc finger protein KOX5, is a 682 amino acid nuclear protein that contains one KRAB domain and eighteen C₂H₂-type zinc fingers. ZNF45 is expressed in many tissues including pancreas, colon small intestine, heart, brain, testis, ovary, prostate, skeletal muscle, placenta, thymus, blood leukocytes and spleen.

REFERENCES

1. Constantinou-Deltas, C.D., et al. 1992. The identification and characterization of KRAB-domain-containing zinc finger proteins. *Genomics* 12: 581-589.
2. Constantinou Deltas, C., et al. 1996. Complete coding sequence, exon/intron arrangement and chromosome location of ZNF45, a KRAB-domain-containing gene. *Cytogenet. Cell Genet.* 75: 230-233.
3. Shannon, M., et al. 1998. Tandem zinc-finger gene families in mammals: insights and unanswered questions. *DNA Seq.* 8: 303-315.
4. Shannon, M. and Stubbs, L. 1998. Analysis of homologous XRCC1-linked zinc-finger gene families in human and mouse: evidence for orthologous genes. *Genomics* 49: 112-121.
5. Tang, M., et al. 2002. Zinc finger gene clusters and tandem gene duplication. *J. Comput. Biol.* 9: 429-446.
6. Hall, T.M. 2005. Multiple modes of RNA recognition by zinc finger proteins. *Curr. Opin. Struct. Biol.* 15: 367-373.
7. Trivedi, A.K., et al. 2007. Proteomic identification of C/EBP-DBD multiprotein complex: JNK1 activates stem cell regulator C/EBP α by inhibiting its ubiquitination. *Oncogene* 26: 1789-1801.
8. Bonetti, A., et al. 2009. A follow-up study of chromosome 19q13 in multiple sclerosis susceptibility. *J. Neuroimmunol.* 208: 119-124.
9. Online Mendelian Inheritance in Man, OMIM™. 2010. Johns Hopkins University, Baltimore, MD. MIM Number: 194554. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: ZNF45 (human) mapping to 19q13.31.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

ZNF45 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ZNF45 shRNA Plasmid (h): sc-97422-SH and ZNF45 shRNA (h) Lentiviral Particles: sc-97422-V as alternate gene silencing products.

For independent verification of ZNF45 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-97422A, sc-97422B and sc-97422C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ZNF45 siRNA (h) is recommended for the inhibition of ZNF45 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ZNF45 gene expression knockdown using RT-PCR Primer: ZNF45 (h)-PR: sc-97422-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.