

GLTSCR1 siRNA (h): sc-97514

BACKGROUND

GLTSCR1 (glioma tumor suppressor candidate region gene 1 protein) is a 1,560 amino acid protein. A polymorphism within the GLTSCR1 gene has been associated with prostate carcinoma aggressiveness. GLTSCR1 is expressed at moderate levels in heart, brain, placenta, skeletal muscle and pancreas, with lower levels in lung, liver and kidney. GLTSCR1 is highly conserved among humans, chimps, mice and rats. GLTSCR1 exists as two alternatively spliced isoforms and is encoded by a gene that maps to human chromosome 19q13.33. Chromosome 19 consists of over 63 million bases, houses approximately 1,400 genes and is recognized for having the greatest gene density of the human chromosomes. It is the genetic home for a number of immunoglobulin (Ig) superfamily members, including the killer cell and leukocyte Ig-like receptors, a number of ICAMs, the CEACAM and PSG families and Fc receptors (FcRs).

REFERENCES

1. Teglund, S., Olsen, A., Khan, W.N., Frängsmyr, L. and Hammarström, S. 1994. The pregnancy-specific glycoprotein (PSG) gene cluster on human chromosome 19: fine structure of the 11 PSG genes and identification of 6 new genes forming a third subgroup within the carcinoembryonic antigen (CEA) family. *Genomics* 23: 669-684.
2. Smith, J.S., Tachibana, I., Pohl, U., Lee, H.K., Thanarajasingam, U., Portier, B.P., Ueki, K., Ramaswamy, S., Billings, S.J., Mohrenweiser, H.W., Louis, D.N. and Jenkins, R.B. 2000. A transcript map of the chromosome 19q-arm glioma tumor suppressor region. *Genomics* 64: 44-50.
3. Online Mendelian Inheritance in Man, OMIM[™]. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 605690. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/605690>
4. Yin, J., Rockenbauer, E., Hedayati, M., Jacobsen, N.R., Vogel, U., Grossman, L., Bolund, L. and Nexø, B.A. 2002. Multiple single nucleotide polymorphisms on human chromosome 19q13.2-3 associate with risk of Basal cell carcinoma. *Cancer Epidemiol. Biomarkers Prev.* 11: 1449-1453
5. Yang, P., Kollmeyer, T.M., Buckner, K., Bamlet, W., Ballman, K.V. and Jenkins, R.B. 2005. Polymorphisms in GLTSCR1 and ERCC2 are associated with the development of oligodendrogliomas. *Cancer* 103: 2363-2372.
6. Barrow, A.D. and Trowsdale, J. 2008. The extended human leukocyte receptor complex: diverse ways of modulating immune responses. *Immunol. Rev.* 224: 98-123.

CHROMOSOMAL LOCATION

Genetic locus: BICRA (human) mapping to 19q13.33.

PRODUCT

GLTSCR1 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GLTSCR1 shRNA Plasmid (h): sc-97514-SH and GLTSCR1 shRNA (h) Lentiviral Particles: sc-97514-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GLTSCR1 siRNA (h) is recommended for the inhibition of GLTSCR1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

GLTSCR1 (H-10): sc-515086 is recommended as a control antibody for monitoring of GLTSCR1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GLTSCR1 gene expression knockdown using RT-PCR Primer: GLTSCR1 (h)-PR: sc-97514-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.